Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

COBRA® Electrosurgical Unit

Catalog Number 604-4810B, 604-4811B, 604-4810BR, & 604-4811BR

TABLE OF CONTENTS

ENGLISH	2
DUTCH	49
FRENCH.....	100
GERMAN.....	150
ITALIAN	203
POLISH.....	254
PORTUGUESE	308
SPANISH.....	360

Models
4810B/4810BR, 4811B/4811BR

Instructions for Use

Manufactured by:
AtriCure Incorporated
7555 Innovation Way
Mason, Ohio 45040 USA
Customer Service:
1-866-349-2342 (toll phone)
1-513-755-4100 (phone)

AtriCure Europe B.V.
De entree 260
1101 EE Amsterdam
The Netherlands
+31 20 7005560
ear@atricure.com

ENGLISH

1. COMPONENT AND ACCESSORIES..... 4
2. INDICATIONS AND CONTRAINDICATIONS 6
3. WARNINGS, PRECAUTIONS, AND ADVERSE REACTIONS..... 6
4. GETTING STARTED..... 10
5. CONTROLS, DISPLAYS, AND ACCESSORIES..... 10
6. ACCESSORIES..... 20
7. OPERATIONAL SEQUENCES 20
8. CONTROLLING THE ESU 22
9. INSTRUCTIONS FOR USE..... 27
10. TROUBLESHOOTING THE ESU 27
11. SERVICE AND MAINTENANCE 29
12. PRODUCT DISPOSAL 31
13. REPAIR AND RETURN POLICY 31
14. TECHNICAL SPECIFICATIONS 42
15. END OF LIFE..... 47
16. LIMITED WARRANTY AND DISCLAIMER 47
17. CONTACT INFORMATION 48

CAUTIONS

Read all warnings, cautions, precautions and instructions provided with this unit before using.

Read the instructions, warnings, and cautions provided with electrosurgical accessories before using. Specific instructions for accessories are not included in these Instructions for Use. When handling this product, also pay close attention to all “Danger,” “Warning,” and “Caution” statements given in these instructions.

SAFETY

The safe and effective use of the COBRA Electrosurgical Unit (ESU) depends to a large degree on factors solely under the control of the operator and medical staff being properly trained. These instructions and equipment are for use only by qualified medical professionals trained in the particular technique and surgical procedure to be performed. It is intended as a guide for using COBRA ESU system. It is important that users read, understand and follow the instructions supplied with the unit.

The surgeon should be familiar with the medical literature, complications and hazards of using RF energy in cardiac procedures.

Any service needs or repairs of the unit must be performed by a trained authorized AtriCure service representative only.

To promote the safe use of the ESU system, this section presents the warnings and cautions that appear throughout these instructions. It is important that you read, understand and follow the instructions in these warnings and cautions so that you can operate this system with maximum safety. It is also important that you read, understand, and follow the instructions for use entirely herein this manual.

1. COMPONENT AND ACCESSORIES

System Description

The COBRA® Surgical System (see Figure 1-1) consists of the COBRA Electrosurgical Unit (ESU), Model 604-4810B, 604-4810BR (115 VAC) Hospital-grade power cord, or 604-4811B, 604-4811BR (230 VAC) Hospital-grade power cord, a COBRA® RF Ablation Probe, and accessories.

ESU

The ESU supplies radiofrequency energy in the range of 450-470 kHz. It operates in a temperature controlled, power limited manner, based on operator settings and temperature feedback provided by thermocouples in the Probe. See Section 5 for more information regarding the ESU and its controls, displays, and accessories (Instrument Cable and ESU Remote Cable).


RF Ablation Probes

Probes are sterile, single use electrosurgical devices intended to be used with the ESU. Probes contain independently selectable electrodes on its distal section. Only Estech COBRA® Probes may be used with the ESU. Refer to the specific probe IFU for further instructions.

Additional Accessories

To avoid incompatibility and unsafe operation, use the following accessories with the ESU system.

- Covidien (Valleylab) Patient Return Electrode, Reference Number E7506
- COBRA Instrument Cable, Model Number 659
- Probe cooling subsystem (optional)

Figure 0-1: ESU System**Legend**

1. **Electrosurgical Probe and Instrument Cable**
2. **Electrosurgical Unit**
3. **Disposable Indifferent (Dispersive) Pad (DIP) Electrodes**

Principles of Operation

Two thermocouples are located at each of the Probe electrodes. These thermocouples measure the temperature at the electrode. The temperature signals from the electrodes are analyzed, and RF current is distributed to selected electrodes to maintain the ESU operator set temperature.

Note: In this manual the Model 604-4810B/604-4810BR/604-4811B/604-4811BR COBRA® ElectroSurgical Unit together with the accessories shown in Figure 1-1 are referred to as the 'System', the Model 604-4810B/604-4810BR/604-4811B/604-4811BR ElectroSurgical Unit alone is referred to as the 'ESU' and the various COBRA® RF Ablation Probes are referred to as the 'Probe'.

2. INDICATIONS AND CONTRAINDICATIONS

Indications for Use

The Estech ESU is intended to supply controlled RF energy to the electrodes of Estech RF Ablation Probes. Please see the specific Probe's Instructions for Use for the Probe's intended use.

Contraindications

- This system is contraindicated for tissue coagulation in any situation where, in the physician's opinion, excessive thermal damage to tissue, or collateral damage to adjacent tissue not intended for coagulation may result.
- Local or systemic infection

3. WARNINGS, PRECAUTIONS, AND ADVERSE REACTIONS

Before operating the System, read the following warnings carefully:

This equipment is for use only by trained, licensed physicians.

The voltage selector must be set to the same setting as the power entry module.

Grounding reliability can only be achieved when the power supply cord is connected to a receptacle marked Hospital Only or Hospital Grade.

Connect the power cord to a properly polarized and grounded power source with the frequency and voltage characteristics that match those on the back panel of the unit.

Always turn off and unplug the unit before cleaning to prevent electric shock hazard.

For patient safety use the unit only if the self-test has been completed as described.

Only use the ESU with Estech ablation probes. See www.atricure.com for available Estech ablation probes.

Do not use Estech ablation probes or cables that are damaged in any way.

Do not use the ESU in the presence of fire and explosion hazards in the operating room:

- Flammable substances
- Flammable gases
- Oxygen enriched atmospheres
- Oxidizing agents

The sparking of heating associated with electrosurgery can provide an ignition source.

Use ESU with caution in the presence of internal and external devices such as pacemakers or pulse generators. Interference produced by the use of electrosurgical devices can cause pacemakers to enter an asynchronous mode or can block pacemaker effect entirely.

If patient has an implantable cardioverter defibrillator (ICD), contact the ICD manufacturer for instructions before performing procedure.

Do not turn unit volume down to an inaudible level. This may prevent the surgical team from being aware of alarm and alert messages.

For RF ablation procedures where the high frequency current could flow thorough parts of the body having a relatively small cross-sectional area, the use of bipolar techniques may be desirable in order to avoid unwanted tissue damage.

Do not position the unit in a way that prevents cord removal. The cord acts as an emergency power cut-off.

The equipotential stud should only be connected to IEC 60601-1 compliant equipment to avoid excessive leakage currents.

PRECAUTIONS:

Observe these precautions before using the System:

Do not attempt to operate the System before thoroughly reading this Operator's Manual. It is important that the equipment's operating instructions be read, understood and followed. For future reference, retain this Manual in a convenient, readily accessible place.

When used with accessories that require fluid cooling, the ESU generator should be set up so that it is not in proximity to the fluid cooling subsystem. **THE ESU IS NOT WATERPROOF.**

The operator should keep the set temperature and power limit as low as possible to achieve the desired end effect. This avoids excessive thermal damage to tissue, or collateral damage to adjacent tissue not intended for coagulation.

Read and follow the Disposable Indifferent (Dispersive) Pad (DIP) Electrode manufacturer's instructions for use. Use only DIP Electrodes which meet or exceed ANSI/AAMI IEC 60601 requirements.

Placement of DIP electrodes on thighs can be associated with higher impedances, which can result in automatic RF power shut-off.

The ESU is capable of delivering significant electrical power. Patient or operator injury can result from improper operation of the System. During power delivery, the patient should not be allowed to come in contact with

metal surfaces which are earthed or have appreciable capacitance to earth.

- Apparent low power output or failure of the equipment to function correctly at normal settings may indicate faulty application of the DIP Electrode or failure of an electrical lead. ***Do not increase settings before checking for obvious defects or misapplication.***

- The risk of igniting flammable gases or other materials is inherent in the application of RF power. Precautions must be taken to restrict flammable materials from the surgery site.
- Electromagnetic interference (EMI) produced by the ESU during the delivery of RF power may adversely affect the performance of other equipment.
- DIP Electrodes should be reliably attached with their entire area to the patient's body and as close to the operating field as possible.
- Skin-to-skin contact (for example, between the arms and the body of the patient) should be avoided, for example by insertion of dry gauze.
- When the System and physiological monitoring equipment are used simultaneously on the same patient, monitoring electrodes should be placed as far as possible from the Electrosurgical Probe. Needle monitoring electrodes are not recommended. In any case, monitoring systems incorporating high frequency current limiting devices are recommended.
- The Instrument Cable to the Electrosurgical Probe should be positioned in such a way that contact with the patient or other leads is avoided.
- Failure of the ESU could result in unintended increase of output power.
- Pacemakers and implantable cardioverter/defibrillators can be adversely affected by RF signals. Refer to the manufacturer's Directions for Use.
- Regularly evaluate and inspect re-usable cables and accessories.
- AtriCure relies on the physician to determine, assess and communicate to the individual patient all foreseeable risks of the electrosurgical procedure.

POTENTIAL ADVERSE REACTIONS:

The following adverse events are possible:

- Excessive thermal damage to tissue
- Collateral damage to adjacent tissue not intended for coagulation.
- Wrapping the accessory cord or return electrode cords around metal objects. This may induce current that could lead to shocks, fires or injury to patient or surgical team.
- The potential for alternate site burns. To avoid this, do one or more of the following:
 - Position the return electrode to provide direct current route between the surgical site and the return electrode which avoids skin-to-skin contact areas.

- Place patient return electrode's according to the manufacturer's instructions.

4. GETTING STARTED

Contents


The ESU's shipping carton contains the components identified below:

- 1ea – ESU
- 1ea– Power Cord
- 1ea – System Operator's Manuals

To unpack the ESU, follow these steps:

1. Unpack the ESU carefully and inspect visually for damage.
2. Notify the carrier immediately if the shipment carton is damaged.
3. Verify that the shipping package contains all items listed above and on the packaging list are received.


If the unit or any accessories are damaged, notify AtriCure's Customer Service immediately. Do not use any damaged equipment. Refer to Section 13 for returning the shipment to AtriCure.

5. CONTROLS, DISPLAYS, AND ACCESSORIES

Front Panel

The operator controls for the System are located on the front panel of the ESU. The buttons on the ESU's front panel allow the user to control the operation of the ESU. The lights/LEDs on the front panel show the ESU's status and/or report error conditions. This section describes the operator controls on the front panel (shown in Figure 5-1).

Figure 5-1: ESU Front Panel**Power Display**

The POWER Display (see Figure 5-2) shows the RF power limit setting or the actual output in Watts. When the ESU is initially powered on and has passed its POST (Power On Self Test) the POWER Display displays a simple pair of bars on the rightmost digits: “- -”.

When a probe is connected, the ESU detects the probe and displays a “0” in the POWER display. When the user selects electrodes for energy delivery, the POWER display will show the Power Limit setpoint.

After the user presses the RF POWER CONTROL Button and energy delivery has begun, the POWER Display will show the actual power that the Probe is delivering into the tissue.

The displayed output power may be lower than setpoint power depending on the measured tissue temperature, the tissue impedance and the number of electrodes that are currently active. When RF power is discontinued, the POWER display reverts to the power limit value.

Figure 5-2: POWER Display and Buttons

POWER SET Buttons


The user can adjust the power limit setting up or down from the default settings via the two POWER SET buttons   located below the POWER display.

Timer Display

The TIMER displays the selected treatment time setpoint and during energy delivery, the remaining treatment time in seconds (see Figure 5-3). When the ESU is initially powered on and has passed its POST (Power On Self Test) the TIMER display displays a bar on each of the digits: “- - -”.

When a probe is connected, the ESU detects the probe and displays the probe’s default Timer setpoint in the TIMER display. While the system is in STANDBY, the operator can adjust the treatment time using the Timer Buttons. Once energy delivery has begun, the TIMER displays the remaining time and it is no longer possible for the user to change the Timer setpoint.

Figure 5-3: TIMER Display and Buttons

Timer SET Buttons

The TIMER set (increase/decrease) Buttons   allow the operator to select the treatment time for the Deliver mode.

TEMPERATURE Display


The TEMPERATURE display (see Figure 5-4) shows the temperature setpoint (desired tissue temperature) in degrees Centigrade. When the ESU is initially powered on and has passed its POST (Power On Self Test) the TEMPERATURE display displays a bar on each of the digits: “- - -”.

When a probe is connected, the ESU detects the probe and displays the probe’s default Temperature setpoint in the TEMPERATURE display. While the system is in STANDBY, the operator can adjust the treatment temperature using the Temperature Buttons. Once energy delivery has begun, the TEMPERATURE Display displays the current electrode temperature and it is no longer possible for the user to change the Timer setpoint.

When the ESU is in the RF DELIVER Mode, the TEMPERATURE Display indicates “LO” if the measured temperature is less than 15° C, “HI” if

temperature is greater than 120° C. Otherwise, it shows the measured temperature.

Figure 5-4: TEMPERATURE Display and Buttons

TEMPERATURE SET Buttons

The TEMPERATURE SET (increase/decrease) Buttons   select the temperature that the ESU attempts to maintain during RF power delivery. During Deliver Mode, the ESU automatically adjusts power, within the user-selected upper power limit, to achieve the desired tissue setpoint temperature.


SET/MEASURED LED Indicators

The SET and MEASURED LEDs (see Figure 5-5) are located to the right of the numeric displays and indicate what is currently being displayed in the POWER, TIMER and TEMPERATURE displays.

During RF DELIVER mode, the yellow MEASURED Indicator is normally lit and the numeric displays indicate measured power, remaining time and current temperature.

Whenever any of the up/down buttons for the POWER, TIMER or TEMPERATURE is pressed, the yellow MEASURED LED turns off. The green SET LED then turns on and the numeric display their respective setpoint values.

Figure 5-5: SET and MEASURED LEDs

RF POWER CONTROL Button/Indicator

The RF POWER CONTROL button (see Figure 5-6) contains an integrated yellow Indicator light. It controls the two operational modes of the ESU:

- STANDBY mode – This mode does not deliver RF Energy to the attached Probe and allows the operator to adjust the setpoint parameters.
- RF DELIVER mode – In this mode the ESU delivers RF power to the attached Probe.

Figure 5-6: RF POWER CONTROL Button and Indicator

When the ESU is initially powered on but no Probe is connected, the RF POWER CONTROL Indicator is not lit.

When the ESU detects that a probe is connected, the RF POWER CONTROL Light flashes, indicating the ESU is in the STANDBY mode and prepared for RF power delivery.

Pressing the yellow RF POWER CONTROL Button once places the ESU into the RF DELIVER mode and allows the probe to transfer RF power to the tissue. In this mode, the RF POWER CONTROL Indicator is lit continuously.

When the treatment is completed, either by the Timer expiring or the operator manually turning off the RF energy, the ESU returns to the STANDBY mode and the RF POWER CONTROL Light resumes blinking.

When the user initially presses the RF POWER CONTROL Button to enter the RF DELIVERY mode, a quick test of temperature sensor functionality is made prior to delivery of RF power. If the ESU detects a defective probe with faulty temperature sensor(s) on more than one electrode the ESU does not enter the DELIVER mode. It re-enters the STANDBY mode and illuminates the REPLACE PROBE LED, indicating a faulty probe has been detected. Also, the ELECTRODE SELECTION LED(s) for the defective electrode(s) blinks.

ELECTRODE SELECTION Buttons

The ESU is capable of delivering RF power and maintaining a desired set temperature on up to seven electrodes simultaneously. Seven ELECTRODE SELECTION Buttons (see Figure 5-7) located on the ESU Front Panel are used to select one to seven electrodes, depending on probe configuration, for simultaneous RF delivery.

The number of available electrodes which may be selected depends on the probe type. At least one electrode must be selected to enable RF delivery; otherwise, the ESU remains in STANDBY and emits the Error tone if the RF ON Button is pressed.

A defective probe electrode may not be selected. If a probe electrode is not selectable, the ESU emits the ERROR tone when the corresponding ELECTRODE SELECTED Button, immediately below the indicator, is depressed and the SELECTED ELECTRODE LED will not illuminate.

Figure 5-7: ELECTRODE SELECTION Buttons and LEDs

Depressing an ELECTRODE SELECTED Button alternately enables or disables that electrode. During delivery of RF DELIVER mode, the ELECTRODE SELECTED Buttons are not operable; they can only be used in the STANDBY Mode.

SELECTED ELECTRODE Indicators

Each of the seven SELECTED ELECTRODE green indicators (see Figure 5-7) displays the status of a probe electrode. If the Indicator is continuously illuminated, the electrode is enabled for RF power delivery. If the Indicator is NOT illuminated, the electrode is disabled and will not deliver RF power.

During the DELIVER Mode, the Selected Electrode Indicators display those electrode(s) being used to deliver RF power.


Electrode(s) in which problems have been detected will have a corresponding flashing Indicator.

ISOLATED PATIENT Connectors

NOTE: The 40-pin ISOLATED PATIENT CONNECTOR is provided for future use only.

The ISOLATED PATIENT CONNECTORS (see Figure 5-8) provide for connection of the Instrument Cable and/or Electrosurgical Probe to the ESU.

Two connectors are available. On the left is a 30 position connector. On the right is a 40 position connector. A sliding door covers the unused connector. The connectors have locking mechanisms for secure cable installation and are keyed for proper alignment of the mating plug.

Figure 5-8: Isolated Patient Connectors (right connector exposed)

The Instrument Cable and Electrosurgical Probe must be installed into one of the ISOLATED PATIENT CONNECTIONS before the ESU can deliver RF power.

CHECK/REPLACE Indicators

The green CHECK and yellow REPLACE Indicators are located beside the two Probe Connectors (

Figure 5-9) and indicate the status of the probe currently connected to the ESU. If no probe is connected or is not detected by the ESU, the CHECK LED is lit. If a defective probe is detected, the REPLACE LED is lit.

Figure 5-9: CHECK and REPLACE Indicators

INDIFFERENT Electrode Connectors and Indicators

The ESU provides two front panel connectors for Indifferent Electrodes. Each connector has an associated green Indicator located directly above the connector (see Figure 5-10).

During RF DELIVER Mode, the ESU monitors the RF current through each electrode. If the measured current on an Indifferent connector rises close to the safe limit, the corresponding indicator will start to blink. If the current rises above the safe threshold, the ESU will detect this condition as an error and terminate the RF DELIVER mode.

Figure 5-10: Indifferent Electrode Connectors and Indicators**SYSTEM FAULT Indicator**

The ESU provides a red SYSTEM FAULT indicator on the front panel that illuminates whenever the internal controller detects a system fault that requires user intervention (see Figure 5-11).

Figure 5-11: SYSTEM FAULT Indicator**Rear Panel**

This section describes the output ports and power controls on the rear panels shown in Figure 5-12 through Figure 5-15.

POWER Switch

The power switch (rocker switch) on the rear panel powers the ESU ON (“1” position) or OFF (“0” position).

Volume Control Knob

This knob provides the user with volume control for the audio signal that accompanies RF power delivery. If no audio tone is heard, check the knob to determine if the volume level has been adjusted too low.

Power Cord

The Power Cord is “Hospital Grade” and must be installed at an AC electrical wall outlet designated “Hospital Grade” or “Hospital Only.”

MAINS Fuse Holder

Housing for AC fuse. (See ESU rear panel for appropriate type of fuse.)

Isolated Serial Port

This port is labeled “RS-232.” It is provided for use with future optional ESTECH approved accessories only.

Ground Stud

The Protective Earth Ground stud provides a chassis ground connection for tying the ESU into a safety grounding point.

Figure 5-12: Model 604-4810B Rear Panel

Figure 5-13: Model 604-4810BR ESU Rear Panel

Figure 5-14: Model 604-4811B ESU Rear Panel

Figure 5-15: Model 604-4811BR ESU Rear Panel

6. ACCESSORIES

Instrument Cable

The Instrument Cable is inserted into one of the ISOLATED PATIENT CONNECTORS (see Figure 5-8) and connects an Electrosurgical Probe to the ESU. The mating plug connectors have a locking mechanism for secure cable installation.

Electrosurgical Probe

The ESU supports monopolar Electrosurgical Probes with any number of electrodes. However, the ESU is only capable of driving up to 7 electrodes simultaneously.

The ESU is also designed to accept monopolar Electrosurgical Probes that have an integral Instrument Cable. Probes with integral cables are designed for single use applications.

Indifferent Electrode

The ESU supports two Indifferent Electrodes that can be connected to connectors on the front panel.

7. OPERATIONAL SEQUENCES

In-service

An initial orientation may be conducted with our AtriCure account representative. All set-up steps and function checks will be covered during your in-service. The features and functionality of your unit will be also covered at that time.

Initial Installation

	 CAUTION
<i>Follow the instructions in this section carefully to prepare the System for operation.</i>	

1. Connect the ESU Power Cord plug into a properly grounded AC electrical outlet designated "Hospital Grade" or "Hospital Only".
2. Connect the ESU Power Cord plug into a properly grounded AC electrical outlet designated "Hospital Grade" or "Hospital Only".

CAUTION: *Never use an outlet without a grounding connection.*

3. Position the ESU for easy viewing of its front panel displays.
4. Install two DIP (Dispersive Indifferent Patch) Electrodes into the INDIFFERENT ELECTRODE receptacles located on the ESU front panel (see Figure 0-1) as follows:

- a. Read the manufacturer's manual before installing the DIP Electrode pads.
- b. Gently push each DIP Electrode fitting straight and firmly into place.
- c. To disconnect a DIP Electrode, grasp the DIP Electrode fitting and gently pull it out from the ESU front panel receptacle.

CAUTION: Do not disconnect a DIP Electrode Connector by pulling on its cable.

- d. Before use, check the DIP Electrodes for damage to the sealed foil packages. Exposure to air, due to a damaged package, could cause the DIP Electrodes to become dry and limit their grounding capability.
- e. Be sure the pads are moist and sticky to the touch before placing on the patient.

NOTE: Do not attempt to relocate the patient grounding pad after initial application. Electrode gel is NOT required and should NOT be used.

- f. Two DIP Electrodes should be placed on a well-vascularized, convex skin surface that is in close proximity to the electrosurgery site. Avoid scar tissue, bony prominence, adipose tissue, or any areas where fluid may pool. Shave, clean, and dry the application site as needed. Avoid wrinkles or folds when applying the pads as wrinkles and folds decrease conductivity.

NOTE: The ESU limits RF output to 1 amp of current flow through either of the DIP Electrode connectors.

5. Move the sliding door to access the appropriate receptacle on the ESU front panel.
6. When the Electrosurgical Probe is ready for use, install the Instrument Cable plug into the appropriate ISOLATED PATIENT CONNECTOR. Plug the cable in by carefully lining up the spline on the connector with the top of the socket. Push the connector firmly into the socket until it locks in place.
7. If using a Probe with a Reusable Instrument Cable, connect the Electrosurgical Probe to the far end of the Instrument Cable.
8. To disconnect the Instrument Cable from the ESU, grasp the connector and pull it straight out of the socket.

CAUTION: Do not disconnect the Instrument Cable by pulling on the cable.

Powering On the ESU

1. Turn the ESU on by pressing the power switch (to the “1” position) located in the power input module on the ESU rear panel (Figure 5-12 through Figure 5-15).
2. The ESU automatically initiates a self-test procedure which is indicated by the illuminated front panel and ESU Remote displays and continuous audio tone generated for approximately 2 seconds. If no system malfunction is detected, the ESU proceeds to the STANDBY Mode.

NOTE: *The Electrosurgical Probe is not required to be connected to the ESU to complete the ESU Power On Self-Test (POST).*

3. If a system malfunction is detected during self-test, all front panel displays remain lit and the ESU does not operate. If any malfunctions are found during self-test, power the ESU OFF, then back ON to attempt to clear the malfunction. The self-test will repeat.

CAUTION

Call AtriCure for Service. A second self-test failure is indicative of an ESU malfunction and should be referred to AtriCure for service. The ESU does NOT operate unless the initial self-test has been successfully completed.

Turning the RF Power Delivery “ON”

To deliver RF power to the Probe, follow these steps:

1. Ensure that all connections are secure and correct.
2. Verify that all preceding selections have been made.

To begin RF power delivery, press the RF POWER CONTROL Button once. The POWER Display shows the RF power delivered to the Probes (in watts).

Turning the RF Power Delivery “OFF”

To discontinue RF power delivery manually during operation (DELIVER Mode), Press the RF POWER CONTROL Button.

8. CONTROLLING THE ESU

ESU in STANDBY Mode

*Refer to **Figure 5-1** for the location of displays and buttons.*

1. When the ESU has successfully completed POST, the ESU enters the STANDBY mode. The RF POWER CONTROL Light will blink to indicate that the ESU is in STANDBY mode.

2. The STANDBY mode is designed for setting the desired parameters for RF power delivery. It is automatically initiated after completion of a successful ESU self-test. STANDBY mode is indicated by a flashing RF POWER CONTROL Button on the ESU front panel; TEMPERATURE, TIME, and POWER Displays are continuously lit.

NOTE: *The numeric displays will show only bars: “- - -” in STANDBY when no probe is connected to the ESU.*

3. Electrodes cannot be selected and the Power, Time and Temperature setpoints cannot be adjusted without first attaching a Probe. Depressing any front panel control button will have no effect until a valid Probe is connected.
4. When the user attaches a Probe, the ESU determines the default setpoints for POWER, TIME and TEMPERATURE for the Probe and those values are displayed on the appropriate front panel displays.
5. Once a Probe is connected, the operator can select electrodes for RF delivery, raise or lower the default Power Limit, Treatment Time or Temperature as required.

ESU in RF DELIVER Mode

1. After the Probe settings have been configured, the user can press the RF POWER CONTROL Button once to place the ESU in the RF DELIVER Mode.
2. When the System is in RF DELIVER Mode, the POWER CONTROL Indicator will remain lit continuously.
3. The POWER display shows (in watts) RF power delivered to the Probe. The value may fluctuate slightly due to changes in tissue contact.
4. The TEMPERATURE display shows (in °C) the maximum temperature measured by any selected Probe electrode. The value may be slightly lower than the actual maximum tissue temperature. The TEMPERATURE Display will indicate “LO” if the measured temperature is less than 15°C, “HI” if temperature is greater than 120°C. Otherwise, it displays the measured temperature.
5. The TIMER Display will indicate the amount of time remaining until the unit automatically shuts down. The Timer will not actually start counting until all of the selected electrodes have reached 50°C. When loss of contact with tissue is detected at all Probe electrodes, the TIMER Display stops for up to 10 seconds. If contact with tissue is not restored within 10 seconds, the ESU returns to STANDBY, emits the Error tone and RF delivery ceases.
6. During RF DELIVER Mode, the TEMPERATURE Display will indicate the current temperature of each of the selected Probe electrodes.
7. RF DELIVER Mode will automatically terminate when the Timer decrements to zero. The operator can manually terminate RF

DELIVER Mode at any time by pressing the RF POWER CONTROL Button once.

About the Power Limit

During RF delivery, the ESU adjusts the power output to maintain the actual measured temperatures of each selected electrode at the temperature setpoint.

The RF Power Limit (150 Watts maximum) is the maximum amount of power that the ESU will deliver to achieve and maintain the desired tissue temperature setpoint. The ESU only supplies as much RF power as necessary to achieve the desired tissue temperature.

If the RF power limit is equal to or greater than the amount of power necessary to achieve and maintain the desired tissue temperature, then the ESU supplies the appropriate amount of power to achieve the temperature setpoint.

If the RF power setpoint is set to a value less than the amount of power required to achieve and maintain the desired tissue temperature, the ESU will not achieve the desired tissue temperature and the measured temperature will remain less than the desired tissue temperature.

Each Probe has a maximum recommended power for each electrode. The default RF Power Limit is the recommended power limit per electrode times the number of selected electrodes. For example, if the probe's per-electrode power limit is 25 Watts (25W), selecting 1 electrode will change the POWER display to "25W", two channels will change the POWER display to "50W", etc.

The ESU allows the operator to adjust the Power Limit, ignoring the default power limit, up or down in 1 Watt increments to a minimum of 1 Watt to a maximum of 150 Watts, regardless of the number of electrodes selected.

Setting the Power Limit

To change the RF power limit (in Watts), follow these steps:

1. Press the appropriate POWER Button   to increase or decrease the RF power limit by 1 Watt.
2. To scroll rapidly to the desired power setpoint, depress and hold the appropriate POWER Button   down.
3. Release the POWER Button   when the POWER Display shows the appropriate RF power limit.
4. During RF DELIVER Mode, the power setpoint can be adjusted in 1-W increments to provide better control during the procedure. While making the adjustment, the SET LED will illuminate to indicate that the value in the POWER Display is the RF power limit setpoint. Once the operator completes the adjustment, the MEASURED LED is then

illuminated to indicate that the value in the POWER Display is once again the current value of measured power.

About the ESU Timer

The ESU has a built-in timer that controls the length of time that RF energy is delivered to the Probe. When a probe is first connected, the TIMER Display will display the default Time of the connected probe.

During STANDBY Mode, the operator can adjust the Treatment time by using the buttons below the time display.

During RF Delivery, the Timer waits until all of the selected electrodes have initially reached 50°C before it starts to decrement down. As it counts down, the TIMER Display shows the remaining Time left for the treatment. When the Timer reaches 0, the ESU discontinues RF energy to the Probe and returns to the STANDBY state.

During RF DELIVER Mode, the Timer cannot be adjusted by the operator.

During RF DELIVER Mode, the operator can manually stop RF energy by pressing the RF POWER CONTROL button once.

Setting the Timer

The default value Timer setpoint for an Electrosurgical Probe connected to the ESU can be adjusted by the operator only during STANDBY mode. During RF DELIVER Mode, the ESU delivers RF Energy to the Probe for the set time period and then returns to STANDBY Mode.

To change the Timer setpoint for the ESU, follow these steps:

1. Ensure that the ESU is in STANDBY mode.
2. Press the appropriate TIMER Button   to increase or decrease the Timer setpoint by 1 second. The SET LED is illuminated to indicate that the value in the TIMER Display is the Timer setpoint.
3. To scroll rapidly to the desired Timer setpoint, depress and hold the appropriate TIMER Button down.
4. Release the TIMER Button   when the TIMER Display shows the appropriate Timer setpoint.

About Temperature Measurement and Control

Each electrode on the Probe has its own temperature sensors. During RF delivery, the ESU monitors those sensors in order to regulate the temperature of each active electrode and the tissue it contacts at the desired temperature setpoint.

During initial heating, the electrodes are significantly below the setpoint and the ESU heats each electrode until it is at the temperature setpoint. The Timer does not start decrementing until all of the selected sensors have initially reached 50°C. Once the sensors reach the set point

temperature, the ESU regulates all of the selected electrodes at the temperature setpoint, turning each off as its temperature rises above the setpoint and turning it back on as it falls below.

For tissue coagulation, it is recommended that the operator initially select the lowest possible set temperature. It should be noted that measured temperature may be slightly lower than maximum tissue temperature. Measured temperature will be influenced by the degree of tissue contact.

Selecting Temperature Setpoint

The default value for desired tissue temperature (temperature setpoint) is determined by the probe type of the Electrosurgical Probe connected to the ESU but can be adjusted by the operator. The ESU then automatically adjusts power, within the user-selected upper power limit, to achieve the desired tissue temperature.

To change the temperature setpoint for the ESU, follow these steps:

1. Ensure that the ESU is in STANDBY mode.
2. Press the appropriate TEMPERATURE Button   to increase or decrease the temperature setpoint by 1°C. The SET LED is illuminated to indicate that the value in the TEMPERATURE Display is the temperature setpoint.
3. To scroll rapidly to the desired temperature setpoint, depress and hold the appropriate TEMPERATURE Button down.
4. Release the TEMPERATURE Button   when the TEMPERATURE Display shows the appropriate temperature setpoint. The MEASURED LED is then illuminated to indicate that the value in the TEMPERATURE Display is the current value of measured temperature.
5. The tissue temperature selection range is 50° to 95°C in increments of 1°C, although the ESU cannot achieve a temperature which is less than the ambient tissue temperature.
6. During RF DELIVER Mode, the temperature setpoint can be increased or decreased by pressing the appropriate TEMPERATURE Button  . When the TEMPERATURE Button   is depressed, the SET LED is illuminated, and the display shows the new temperature setpoint. When the button is released, the MEASURED LED is illuminated, and the display shows the present value of tissue temperature.

About the Audio Control

The user may adjust the volume of the audio signal when the ESU is initially powered on or during a surgical procedure.

Adjusting Audio Volume

1. Adjust the volume of the audio signal by rotating the Volume Control Knob on the ESU rear panel (see Figure 5-12 through Figure 5-15).
2. If desired during a clinical procedure, the Volume Control Knob may be rotated fully counter-clockwise to its lowest setting.

NOTE: *At the minimum setting, the tone may NOT be audible during RF power delivery.*

9. INSTRUCTIONS FOR USE

Once the Electrosurgical Probe has been positioned, RF power is delivered via the ESU resulting in the coagulation of the targeted tissue. Temperature can be monitored during the coagulation procedure through the sensors embedded within Probe electrodes.

1. Follow steps given in the sections, *Initial Installation* and *Powering ON the ESU*.
2. Follow steps given in the *Instructions for Use of the Probe*.
3. Connect the Probe to the ESU by plugging in the Instrument Cable.
4. Select the desired Probe electrodes with the ELECTRODE SELECTED Buttons on either the front panel or the ESU Remote Cable.
5. When the coagulation site has been accessed and the selected Probe Electrodes are in contact against the tissue surface, RF power may be delivered as described in Section 8, *ESU in RF DELIVER Mode*.
6. RF power is delivered to the tissue via the selected probe electrodes, resulting in coagulation under the selected Probe electrodes. Temperature sensors in the electrodes allow temperature control of RF delivery, resulting in more consistent coagulation without over temperature effects.

10. TROUBLESHOOTING THE ESU

Indications of ESU Problems/System Faults


RF power delivery can cease due to operating conditions other than user intervention. To assist the operator, the ESU front panel indicators can provide indications of anomalous operating conditions and equipment setup failures.

Front Panel Fault Indicators:

SYSTEM FAULT Indicator is lit when a non-recoverable software or hardware failure has been detected and the ESU will remain

inoperative. This fault can be triggered by an internal system failure or by operating the ESU outside its safe operating conditions (ambient temperature too high or too low, etc.).

Table 10-1 lists the ESU Indicators and related fault conditions.

Table 10-1: Error Indicators	
Indicator	Description
 CHECK	<ul style="list-style-type: none"> No Probe Connected Invalid Probe Connected Defective Probe Connected Indifferent Electrode Problem Electrode Over-current Condition
 CHECK (blinking)	<ul style="list-style-type: none"> High Load Impedance detected (above 300 ohms)
 REPLACE	<ul style="list-style-type: none"> Probe electrode with 2 defective thermocouples Invalid Probe ID Probe activation has expired Defective Probe
 INDIFFERENT ELECTRODES (GROUNDING PADS) (blinking)	<ul style="list-style-type: none"> Indifferent Electrode in poor contact with patient or tissue (current into the opposite DIP has exceeded 1.1 Amps)
 SYSTEM FAULT	<ul style="list-style-type: none"> Non-recoverable Software or hardware error System operation error

Possible Causes of RF Power Delivery Interruption

RF power delivery can cease due to operating conditions other than user intervention. RF power automatically ceases when:

- The ESU detects a current of 1.1 Amps or more in either of the two DIP Electrode connectors.
- Impedance is less than 25 Ω .
- Impedance is greater than 300 Ω for more than 10 seconds.
- Temperature of a non-selected electrode exceeds 70° C for more than 1 second.
- The ESU detects a current of 0.9 Amps or more in any Probe electrode for more than 3 seconds.
- Temperature of an electrode exceeds 120° C for more than 1 second.
- Temperature of the hottest electrode exceeds the temperature of the coolest electrode by more than 15° C for more than 5 seconds.

System operation errors are detected.

Correcting a CHECK Condition

If the CHECK Indicator LED is lit during RF power delivery, follow these steps:

1. Ensure that the connection between the Probe and ESU is correct and secure.
2. Ensure that the connection between each DIP Electrode and the ESU is correct and secure.
3. Ensure that each DIP Electrode is properly applied.
4. If an improper cable connection is not the cause of the CHECK PROBE condition, then a build-up of coagulum on the Probe electrodes may be the cause, and RF power output should be discontinued.
5. To discontinue RF power delivery, refer to the instructions in Section 7, *Turning the RF Power Delivery "OFF"*.

Correcting a REPLACE Condition

If the REPLACE PROBE LED is lit, follow these steps:

1. Ensure that the connection between the Probe and the ESU is correct and secure.
2. Replace the instrument cable, if the REPLACE PROBE condition persists.
3. Replace the Probe, if replacing the instrument cable does not correct the REPLACE PROBE condition.

Correcting a SYSTEM FAULT Condition

To correct a SYSTEM FAULT Condition:

- 1) Correct any operating conditions that are outside the ESU's normal operating range.
- 2) Turn the ESU "OFF." Wait approximately 4 seconds before turning the ESU "ON". Allow the ESU to perform the Power-On-Self-Test (POST).
- 3) If the SYSTEM FAULT Condition reoccurs, **Call AtriCure for Service**

11. SERVICE AND MAINTENANCE

Servicing the ESU System

The ESU requires no routine service or maintenance. If the ESU fails to operate when plugged into a proper AC power receptacle and the power switch is turned on, check the fuse (refer to Replacing the Fuses in this

chapter).

The ESU contains no user-serviceable parts; disassembly and attempted repair by unqualified personnel may create a hazardous condition and voids the warranty.

Cleaning/Disinfecting the ESU System

NOTE: Do not spray or pour liquids directly on the unit.

NOTE: The unit and/or accessories cannot be sterilized.

WARNING: Ensure Isopropyl Alcohol (IPA) is completely dry before operating the unit.


CAUTION: Avoid caustic or abrasive cleaners

Guidelines

The following guidelines are recommended for cleaning the unit. It is the user's responsibility to qualify any deviations from these processing methods.

1. Disconnect the unit or cart from the outlet before cleaning.
2. If the unit and/or accessories are contaminated with blood or other body fluids, they shall be cleaned before the contamination can dry (within two hours of contamination).
3. The outer surfaces of the unit and/or accessories shall be cleaned with 70% - 90% Isopropyl alcohol (IPA) wipes for a minimum of two minutes. Do not allow fluids to enter the chassis.
4. Pay attention to all areas where fluids or soil may gather, such as under/ around the handles or any tight crevices/ grooves.
5. Dry the unit and/or accessories with a dry, white lint-free cloth.
6. Conduct a final confirmation of the cleaning process by visually inspecting the white cloth for remaining soil.
7. If soil remains on the white cloth, repeat steps 3 through 6.
8. After cleaning is complete, turn the unit on to perform Power On Self-Test (POST). If any errors are received, contact AtriCure to begin return process.

Replacing the Main Fuses

1. Before replacing a fuse in the ESU, disconnect the mains power cord

from the ESU.

2. Pull the fuse holder out of the power entry module. A slotted screwdriver may be used to assist in removing the fuse holder.

Replace the fuse with another of the same type and rating. Refer to the fuse label on the rear panel (see Figure 5-12 through Figure 5-15).

3. Install the new fuse in the fuse holder.

Reinsert the fuse holder in the power input module taking care to install it in the correct orientation.

12. PRODUCT DISPOSAL

After use, dispose of product and packaging in accordance with hospital, administrative and/or local, state, federal and international laws and regulations. Please see section 13 for AtriCure return policy.

13. REPAIR AND RETURN POLICY

AtriCure is responsible for the safety, reliability and performance of the ESU only under the following circumstances:

- Persons authorized by AtriCure has performed assembly, operation, readjustments, modifications or repairs.
- Equipment use is only in accordance with AtriCure Instructions for use.
- The electrical installation of the room complies with the local and regulatory requirements, such as IEC, etc.
- Equipment is to be disposed of or recycled in accordance with local laws.

Please note that infected medical devices must be decontaminated or properly disposed of as medical/biohazard waste.

Before returning the unit, call your AtriCure account representative for assistance. If instructed to return the unit to AtriCure, please obtain a Return Material Authorization number from AtriCure Customer Service. Then clean the Unit and package securely to ensure proper protection of unit. Please record the Returned Material Authorization number on the outside of the container.

Repackaging the ESU

When the ESU requires service or transfer to other location, use the original shipping carton and packing materials to repack and ship the ESU. For shipping, disconnect all accessories, and place the ESU and power cord in the locations reserved for the items in the carton.

14. TECHNICAL SPECIFICATIONS

PERFORMANCE CHARACTERISTICS

Table 14-1 lists the specifications for the System (ESU and Accessories).

Table 14-1: System Specifications	
Description	Specification
Input Power	
Supply Mains/Line Voltage Rating	100 – 120 V~ (Model 4810B) 100 – 120 V~ (Model 4810BR) 220 – 240 V~ (Model 4811B) 220 – 240 V~ (Model 4811BR)
Line Current Rating	4.0 Amps, 50/60 Hz (100 – 120 V~) 2.0 Amps, 50/60 Hz (220 – 240 V~)
Fuse Rating	5.0 A, 50/60 Hz (100 – 120 V~) 2.5 A, 50/60 Hz (220 – 240 V~)
Power Cord Length	10 ft. (3 m)
Instrument Cable	
Length	10 ft. (3 m)
Connector	30-pin Polarized Quick-Connect Connector
Indifferent Electrode Cable	
Current Limit	1.1 Amperes
ESU	
RF Power Outputs	<ul style="list-style-type: none"> RF energy is delivered at 460 kHz, as a sinusoid. 150 W maximum into a minimum impedance of 37.5 Ω, with a 2 A_{RMS} current limit Maximum output voltage is 167

	V_{RMS} . <ul style="list-style-type: none">• Maximum RF power output is based on Probe type.
Impedance	<ul style="list-style-type: none">• Measures 25 Ω and above, if accessing up to 150 W of RF power.• Measures up to 999 Ω.

Table 14-1: System Specifications	
Description	Specification
Temperature	Allows selection of desired tissue temperature (temperature setpoint) within the range of 5° to 95° C.
Normal Operating Condition: With ESU in the DELIVER Mode	ESU measures 15° to 120° C.
Outside of Normal Operating Range	ESU displays “LO” or “HI”.
Power Limit	Allows selection of power limit within the range of 1 W – 150 W.
Time Measurement	0 to 600 seconds in increments of 1 second.
Shutdown Limits	
Measured Temperature	Below 0° C or above 120° C.
Measured Impedance	Below 25 Ω or above 300 Ω for more than 10 seconds.
Measured Current	RF output to the patient is limited if more than 1 A (+20%, -0%) flows in either INDIFFERENT ELECTRODE.
Measurement Accuracy	
Power ¹	$\pm 5\%$, ± 2 W (see Figure 14-1)
Impedance ¹	$\pm 5\%$, ± 5 Ω
Temperature	$\pm 3\%$, $\pm 2^\circ$ C (see Figure 14-2)
Dimensions and Weight	
Height	6 in. (15 cm)
Width	14 in. (36 cm)
Depth	14 in. (36 cm)
Weight	22 lb. (10 kg)
Low-Frequency Leakage Current (50/60 Hz): ²	
Normal Polarity, Intact Chassis Ground	<10 μ A
Normal Polarity, Ground Open	<50 μ A
Reverse Polarity, Ground Open	<50 μ A
Sink Current, at 120 V Applied, All Inputs	<500 μ A
Chassis Source Current, Ground Open	

Notes:

1. Power and impedance measurements are specified at the ESU output connections and do not account for effects that the cable or the probe may have.
2. Source Current, Patient Leads, All Outputs Tied Together

ENVIRONMENTAL SPECIFICATIONS

Table 14-2 lists the environmental specifications for the System.

Table 14-2: Environmental Specifications	
Description	Specification
Storage/Transport Temperature	-40°C to 70°C (-40°F to 158°F) The unit should be gradually returned to the operating temperature range before use and stabilized for one hour before operation.
Relative Humidity	10% to 90%, non-condensing The unit should be gradually returned to the operating environmental conditions before use and stabilized for one hour before operation.
Operation Temperature	10°C to 40°C (50°F to 104°F)
Relative humidity	30% to 75%, non-condensing
Altitude	0 to 2000 meters

LABELING SYMBOLS

Table 14-3: Labeling Symbols	
Image	Description
	Caution
	Dangerous Voltage
	Warning, Electricity
	Defibrillator-Proof Type CF Equipment
	Disposable Indifferent Pad Electrode



Table 14-3: Labeling Symbols	
Image	Description
	Equipotentiality
	Footswitch
	Fuse
	Increase / Decrease
	Loudspeaker Volume
	Mains Port
	Non-ionizing Radiation
	Probe/Instrument with Integrated Cable
	Probe/Instrument Cable
	STOP (RF OFF)
	START (RF ON)
	Select
S232 010101	Serial Data
	Temperature
	Time
	Isolated Patient Circuit

Table 14-3: Labeling Symbols	
Image	Description
	WEEE Label: Contact AtriCure for return and recycling information at the end of useful life
	Non-Sterile
	Manufacturer/Company Address
	Authorized Representative in the European Community
	Product complies with the requirements of directive 93/42/EEC
	Serial Number
	Refer to instruction manual/booklet Note: On equipment "Follow instructions for Use"
	General mandatory action sign
	Danger indicates a hazard with a high level of risk which, if not avoided, will result in death or serious injury.
	Warning indicates a hazard with a medium level of risk which, if not avoided, could result in death or serious injury.
	Caution indicates a hazard with a low level of risk which, if not avoided, could result in minor or moderate injury.
	Catalogue Number
	Quantity
	Date of Manufacture
Rx ONLY	Caution: U.S. Federal law restricts this device to sale by or on the order of a physician or other licensed practitioner.



	<p>Product is certified for both the U.S. and Canadian markets, to the applicable U.S. and Canadian standards.</p>
	<p>Storage temperature range</p>
	<p>Storage humidity range</p>

Table 14-4 IEC Electrical Safety Specifications
Device Description: Class I Equipment, Type CF Patient Applied Parts (Defibrillation-proof)
IP Rating/ Moisture Protection Rating: IPX0 (ESU)
Mode of Operation: Continuous
Electrical Isolation: <ul style="list-style-type: none"> • Leakage current conforms to IEC 60601-1 Dielectric withstanding voltage conforms to IEC 60601-1
EMC Emissions and Susceptibility: The Estech COBRA™ Electrosurgical HF Generator has been tested and found to comply with the limits for medical devices to the IEC 60601-1-2:2005. These limits are designed to provide reasonable protection against harmful interference in a typical medical installation. This system generates, uses, and can radiate radiofrequency energy and, if not installed and used in accordance with the instruction given below, may cause harmful interference to other devices in the vicinity. However, there is no guarantee that interference will not occur in a particular installation.

ELECTRO-MAGNETIC COMPATIBILITY (EMC) SPECIFICATIONS

Table 14-5 IEC EMC Specifications (Emissions)

Guidance and Manufacturer's Declaration – Electromagnetic Emissions		
<p>The Estech COBRA™ Electrosurgical HF Generator is intended for use in the electromagnetic environment specified below. The customer or the user of the Estech COBRA™ Electrosurgical HF Generator should assure that it is used in such an environment.</p>		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF emissions CISPR 11	Group 1	The Estech COBRA™ Electrosurgical HF Generator must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11	Class A	The Estech COBRA™ Electrosurgical HF Generator is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

IMMUNITY EMC SPECIFICATIONS

Table 14-6 IEC EMC Specifications (Immunity)			
Guidance and Manufacturer's Declaration – Electromagnetic Immunity			
The Estech COBRA™ Electrosurgical HF Generator is intended for use in the electromagnetic environment specified below. The customer or the user of the Estech COBRA™ Electrosurgical HF Generator should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 seconds	<5% UT (>95% dip in UT) for 0.5 cycles 40% UT (60% dip in UT) for 5 cycles 70% UT (30% dip in UT) for 25 cycles <5% UT (>95% dip in UT) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Estech COBRA™ Electrosurgical HF Generator requires continued operation during power mains interruptions, it is recommended that the Estech COBRA™ Electrosurgical HF Generator be powered from an uninterruptible power supply or battery.

NOTE <i>UT</i> is the a.c. mains voltage prior to application of the test level.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of typical location in a typical commercial or hospital environment.
Conducted RF IEC 61000-4-6 Radiated RF IEC 61000-4-3	3 Vrms 150 kHz to 80 MHz 3 V/m 80 MHz to 2.5 GHz	3 Vrms 3 V/m	<p>Portable and mobile RF communications equipment should be used no closer to any part of the Estech COBRA™ Electrosurgical HF Generator including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.</p> <p>Recommended separation distance:</p> $d = [1.17] \sqrt{P} \text{ MHz to 800 MHz}$ $d = [1.17] \sqrt{P} \text{ MHz to 800 MHz}$ $d = [2.33] \sqrt{P} \text{ MHz to 2.5 GHz}$ <p>where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range.</p> <p>Interference may occur in</p>

			<p>the vicinity of equipment marked with the following symbol:</p> 
--	--	--	---

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the Estech COBRA™ Electrosurgical HF Generator or any of its components are used exceeds the applicable RF compliance level above, the Estech COBRA™ Electrosurgical HF Generator should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating components or the entire Estech COBRA™ Electrosurgical HF Generator.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Table 14-7 IEC EMC Specifications (Immunity)**Recommended separation distances between portable and mobile RF communications equipment and the Estech COBRA™ Electrosurgical HF Generator**

1. The Estech COBRA™ Electrosurgical HF Generator is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the Estech COBRA™ Electrosurgical HF Generator can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the Estech COBRA™ Electrosurgical HF Generator as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter M		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$	$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$	$d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
50	8.25	8.25	16.50

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

PERFORMANCE GRAPHS

Figure 14-1: Power measurement accuracy graph

Region between the lines indicates possible output range.


Figure 14-2: Temperature measurement graph

Region between the lines indicates possible error range.


Figure 14-3: Power Output versus load at full and half power

Curves start at 25 Ω the minimum impedance required to output RF energy.
 Curves end at 300 Ω the maximum impedance allowed to continuously output RF energy.



Figure 14-4: Peak Output Voltage

15. END OF LIFE

The ESU service life is a maximum 10 years. If the unit performance fails while still under warranty, please contact your local AtriCure sales representative or customer service for a RGA number to return the unit.

16. LIMITED WARRANTY AND DISCLAIMER

Limited Warranties

AtriCure warrants that this product is free from defects in original workmanship and materials. If this product is proved to be defective in original workmanship or original materials, AtriCure, in its absolute and sole discretion, will replace or repair it, less charges for transportation and labor costs incidental to inspection, removal or restocking of product. This limited warranty applies only to original factory delivered products which have been used for their normal and intended uses. AtriCure's limited warranty shall NOT apply to AtriCure products which were installed or calibrated by persons not authorized by AtriCure or which have been repaired, altered, or modified in any way and shall NOT apply to AtriCure products which have been improperly stored or improperly installed, operated or maintained contrary to AtriCure instructions.

Disclaimer and Exclusion of Other Warranties

There are no warranties of any kind which extend beyond the description of the warranties above. AtriCure disclaims and excludes all warranties, whether expressed or implied, of merchantability of fitness for a particular use or purpose.

Limitation of Liability for Damages

In any claim or lawsuit for damages arising from alleged breach of warranty, breach of contract, negligence, product liability or any other legal or equitable theory, the buyer specifically agrees that AtriCure shall not be liable for damages for loss of profits or revenues, loss of use of the product, loss of facilities or services, any downtime costs, or for claims of buyer's customers for any such damages. AtriCure sole liability for damages shall be limited to the cost to buyer of the specified goods sold by AtriCure to buyer which give rise to the claim for liability. The buyer's use of this product shall be deemed acceptance of the terms and conditions of these limited warranties, exclusions, disclaimers and limitations of liability for money damages.

17. CONTACT INFORMATION

Manufactured by:
AtriCure, Incorporated
7555 Innovation Way
Mason, Ohio 45040 USA
Customer Service:
1-866-349-2342 (toll free)
1-513-755-4100 (phone)

European Representative:
AtriCure Europe B.V.
De entree 260
1101 EE Amsterdam
The Netherlands
+31 20 7005560
ear@atricure.com

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Atracure Europe B.V.