# AtriCure®

## frigitronics®


### Operation and Maintenance Manual


# CCS-200

Dual Probe, Non-Electric Cardiac Cryosurgical System

# The ACC2 Cardiac Cryosurgical System

**Front View**

**Rear View**

Nitrous Oxide Cylinders Not Shown

- |  |  |                                      |
|--|--|--------------------------------------|
| 1. Probe Socket - Exhaust (2)          | 9. Battery Compartment (2)               | 17. Filter Assembly - Pressure Gauge |
| 2. Probe Socket - Supply (2)           | 10. Digital Tissue Temperature Indicator | 18. Filter Assembly - Canister       |
| 3. Digital TTI Socket (2)              | 11. Exhaust Pressure Gauge (2)           | 19. Cylinder Valve (2)               |
| 4. Digital Probe Temperature Meter (2) | 12. Digital Temperature Indicator Switch | 20. Cylinder Yoke Assembly (2)       |
| 5. Line Pressure Gauge (2)             | 13. Maximum Freeze Control Valve (2)     | 21. Cylinder Bracket (2)             |
| 6. Digital Timer (2)                   | 14. ON/OFF Control Valve (2)             | 22. Cylinder Strap (2)               |
| 7. Temperature Control Valve (2)       | 15. Main Meter Power Switch              | 23. Mobile Carrier - Casters (4)     |
| 8. Freeze/Defrost Control Valve (2)    | 16. Exhaust Nipple                       |                                      |

## Specifications

### ACC2 System

47" High (119.4cm), 19" Wide (48.25cm)

24" Deep (70cm)

Weight: 74lbs (35.8kg) without cylinders

Note: To convert psi to kg/sq cm, multiply by 0.07.

Two 20lb. Nitrous Oxide Cylinders required

Digital Meters use Four 9Volt Alkaline Batteries

**CAUTION: U.S. Federal law restricts this device to sale by or on the order of a physician.**

## Introduction

The **AtriCure Model ACC2** is a *dual probe*, independently controlled, fully-exhausted, Non-AC Powered (1), self-contained Cardiac Cryosurgical System. This system uses readily available medical grade nitrous oxide gas.

This state-of-the-art instrument was primarily designed for treatment of **cardiac arrhythmias** and can be used for other procedures requiring controlled temperatures down to  $-50^{\circ}$  to  $-60^{\circ}$  depending on the probe used. The system was designed to provide the surgeon with precise temperature control enabling him to perform reversible freezing (2).

The system operates on the Joule-Thomson principal where pressurized gas expands through a fine orifice producing a rapid drop in temperature.

**The system provides the surgeon the capability to:**

- Select either maximum or preset controlled freeze temperatures
- Quickly defrost the probe(s)
- Monitor the probe tip temperature
- Track actual tissue temperature using thermocouple needles
- Time any aspect of the operation, such as the *freeze duration*
- Use two probes consecutively or simultaneously. They can be identical or of different sizes and shapes
- Select from a wide range of probe configurations

## Unpacking and Set-Up

### Unpacking

The ACC2 System includes the following:

- Console with integrated mobile carrier

**The console is supplied with:**

- Four 9 volt batteries
- Replacement “O” rings
- Filter and gauge assembly w/cartridge
- Filter bracket and wrench

**The carrier is supplied with:**

- Four locking casters (23)

Refer to your packing list then check all packing material for accessories and components. Please read this operation manual carefully before using the system.

### Set-Up

1. Install filter bracket and filter assembly (17,18) per the instruction sheet supplied.
2. Install filter cartridge (see To Install Filter Cartridge, Page 7)
3. Install two 9 volt batteries, for the digital display in each of the two battery compartments (9) (Four total). Confirm proper orientation of + and – on 9V battery inside holder. Improper orientation of battery will damage temperature meters.

1. System does not require any external AC power source. Batteries used for Temperature and Timer Display only.
2. Reversible freezing—to cause temporary blockage of the signal when frozen—does not cause tissue necrosis or permanent damage—effect is reversed upon thawing.

## Installation of Cylinders

**Medical grade nitrous oxide** cylinders are locally available worldwide. Consult your yellow pages under Gas Medical or Welding Supplies. Be sure to specify **MEDICAL GRADE** nitrous oxide gas. The ACC2 System is designed for nitrous oxide gas only.

1. Place a full nitrous oxide cylinder against bracket (21) and secure in place by fastening buckle and tightening the straps (22). Repeat procedure for second cylinder. Remove end of cylinder tags to show cylinders are now in use.

### **Note:**

Cylinders must be at room temperature for maximum pressure and optimum performance. Nitrous Oxide at room temperature is approximately 740 psi.

2. Place cylinder yoke nipples (20) into cylinder valves. Hand tighten yoke nut.

**DO NOT OPEN CYLINDER AT THIS TIME.**

3. Two reserve cylinders are recommended. They can be purchased from your local gas supplier.

## Checking The System

### **Gas Exhaust**

The National Institute for Occupational Safety and Health (NIOSH) current recommended exposure limit for nitrous oxide (N<sub>2</sub>O) gas is 25 parts per million (ppm) averaged over the period of administration.

1. Connect the exhaust hose supplied by AtriCure to the console exhaust nipple (16) located on the bottom of the console.
2. Safely vent the other end of the hose in one of the following manners:
  - a. Directly to the outside
  - b. Directly to a **Non-Recirculated** room exhaust ventilation system.
  - c. Directly to a dedicated system for venting anesthetic gases. The flow capacity of the system shall not be less than 40 liters/min.

### **Caution:**

1. Only use the exhaust hose provided by AtriCure Inc.. No substitutes allowed.
2. Do not restrict or kink the exhaust hose in any way.

### **Leaks**

1. Turn *both* console **On/Off** knobs (14) "off".
2. Turn one cylinder (19) "on" by turning one cylinder valve counterclockwise (Gas Gauge (17) should show pressure). Check all exterior connections for leaks.

# Checking the System

## Installation of Probes

- **Do not** remove or install probes unless the line (5) and exhaust pressure gauges (11) read “0” psi.
  - **Do not** restrict, kink, bend, clamp or otherwise damage or restrict probe lines.
  - **Do not** block exhaust holes at the end of probe lines where the probe plugs into the console.
  - **Do not** autoclave or submerge probe in solutions. If autoclaved—*Do Not Use, Return to AtriCure for servicing.*
  - **Do not** attempt to insert or remove probe from console when it is pressurized.
1. Plug desired probe into respective probe sockets (1,2) on front of console. To engage exhaust plug (1), slide retainer ring back on quick disconnect while inserting plug, then release. To engage input plug (2), depress latch while inserting then release. Tug at probe connections to ensure both connectors are latched. Repeat procedure for second probe if required.
  2. Turn the **Main Meter Switch** (15) “on”. The digital probe temperature meters (4) will now display the probe tip temperature.

## **Digital Probe**

### Temperature Meters

The temperature meters (4) are calibrated in degrees Celsius. There is no calibration required. To check if the system is functioning correctly, place the tip (with pressure gauges reading zero) in ice water. The tip temperature display should read approximately 0°C. The temperature meter is for reference only.

### ***Battery Information***

The ACC2 is equipped with three (3) Digital Temperature Meters. These meters are powered by the four 9 volt batteries located in the battery compartment (9) at the rear of the console. Two of the batteries (one each) power the probe meters of System I and System II. The battery clips are marked “Sys I” and “Sys II” respectively.

## Operational Summary

The operation of System I and System II can be done individually or simultaneously. Switch controls slowly to minimize gas surge. The sequence of operation can be changed as required. A flow of gas can be heard in the probe during the operation—this is normal in a functioning probe.

### Operational Set-up

1. With the two probes plugged into the front of the ACC2 console, open *both* cylinder valves (19) fully. Check line pressure gauge (5) to insure adequate line pressure.

**CAUTION: It is not recommended to use the system when the pressure is below 650 psi.**

Use of the system when in the yellow zone will result in reduced performance.

**WARNING: DO NOT OPERATE SYSTEM IN RED ZONE.**

2. Turn the **Freeze/Defrost Control** knobs (8) to the “defrost” position.
3. Close the **Temperature Control** valves (7) by turning fully clockwise. Do not overtighten.
4. Turn the **Maximum Freeze Control** knobs (13) to the “temperature control” position.
5. Turn the **On/Off** control knob (14) to the “on” position.

***You are now ready to start freezing.***

### For controlled Temperature Freezing

- a. Turn the **Maximum Freeze Control** knob to “temperature control”.
- b. *Slowly* turn the **Freeze/Defrost Control** knob to the “freeze position”.
- c. Turn the **Temperature Control** knob counterclockwise to decrease temperature (colder), clockwise to increase temperature (warmer)—monitor the probe temperature meter until the desired tip temperature is obtained. It may be necessary to wait momentarily for the temperature to stabilize before making final adjustments.

### For Maximum Freeze

*Slowly* turn the **Maximum Freeze Control** knob to the “maximum freeze” position.

### To Defrost

The probe can be defrosted or warmed above 0° C at any time. *Slowly* turn the **red Freeze/Defrost Control** knob to the “defrost” position.

## Timers

Your system is equipped with two manual digital timers (6). Timers are constantly “on”. They operate on Lithium batteries, with a life expectancy of 5 years. They may be used to time any portion of the operation, i.e. freeze duration. They are capable of timing up, *elapsed time*, or timing down from a *preset limit*. The timer reads in minutes and seconds. The first two digits are minutes and the second two digits are seconds. The timer can time up to 60 minutes. Read the following instructions carefully then practice using it a few times.

### To Use Your Timer

- A. Minute Mark
- B. Second Mark
- C. LCD Digital Display
- D. Number Keys
- E. Start Button
- F. Stop/Reset Button

### Preset Timing

#### *To Set Time for Countdown*

1. Press the **Stop/Reset** button to set Timer to zero.
2. Press number keys **0** through **9** to set desired minutes and seconds. *Example: To set 12 minutes 38 seconds, press 1,2,3,8.* The first and second keys pressed are the minutes. The third and fourth keys are the seconds.
 

**NOTE:** If you are timing in seconds only, press just two number keys **0** through **9**. *Example: To set 15 seconds, press 1,5.*
3. If the wrong time has been set, simply press the Stop/Reset button to reset time to zero and begin again.
4. When accurate time is set press **Start** button. The seconds will immediately begin counting down.

### When Timing Period is Completed

When the Timer counts down to 00 minutes 00 seconds, an electronic tone will sound. This tone will continue to sound for one minute or until Timer is turned off. (The tone automatically turns off after one minute to extend battery life.) To turn off the tone, simply press any one of the buttons on the Timer: **Start**, **Stop/Reset**, or **0** through **9**.

### Elapsed Timing

#### *To Set Timer for Count-Up*

1. Press the **Stop/Reset** button to set Timer to zero.
2. Press **Start** button once. Timer will immediately begin to count up at 1 second intervals. It will count up to 60 minutes and then automatically recycle itself to zero and start counting up again.
3. To stop the count-up, simply press **Stop/Reset** twice.

### To Stop the Timer During Countdown or Count-Up

If it is necessary to stop Timer while it is counting down or up, press **Stop/Reset** button once and Timer will stop. To start countdown or count-up where you left off, press the **Start** button once. Or, if you want to clear Timer to zero, press **Stop/Reset** button a second time.

To replace the battery in your electronic Timer/Stopwatch, first remove filter assembly from console, then remove bottom cover. Carefully remove Timer backup plate and remove Timer from console. Use a standard Phillips screwdriver to turn the screw on the Timer backside in the direction of the arrow.

Remove the Timer back cover and exhausted battery. Insert a new G-13 size button cell battery making sure the positive (+) side is facing up and replace the Timer back cover. Reverse procedure for installation of Timer.

## Exchanging Components

### To Exchange Probes

1. Turn the On/Off control knobs (14) of System I and/or II “off”, depending on which probe is to be changed.
2. **Do not** remove probe(s) until the system line and exhaust pressure gauge needles touch the silver stop.
3. Remove the probe and exchange. Tug at probe connections to ensure both probe connectors are latched.

### To Exchange Cylinders

1. Turn both **Cylinder Valves** (19) “off” (clockwise). Turn the system **On/Off** knobs “on” and “off” repeatedly. This will relieve the pressure in the system while safely venting the gas.
2. When the cylinder pressure gauge (17) on the filter assembly reads “0” psi, remove cylinder yokes.
3. Exchange Cylinders—Place **On/Off** valves (14) in the “off” position.
4. It is recommended that the filter cartridge be replaced with each change of cylinders. This will help insure clean, trouble-free performance.

**WARNING: USE MEDICAL GRADE NITROUS OXIDE GAS ONLY.**

### To Exchange Filter Cartridge

1. Turn both **Cylinder Valves** (19) “off” (clockwise). Turn the system **On/Off** knobs “on” and “off” repeatedly. This will relieve the pressure in the system while safely venting the gas.
2. When the cylinder pressure gauge reads “0” psi, unscrew the filter canister (18) and remove.
3. Unscrew the filter cartridge and replace with a new cartridge.
4. Inspect filter “O” ring for swelling, nicks, cuts, abrasions, or other physical damage. Replace as necessary.
5. Replace filter canister and tighten securely.

**CAUTION: There should be no gap between canister and housing when tightened.**

6. Turn one **Cylinder Valve** “on” (counterclockwise). Inspect for leaks.

### To Shut Down

1. Turn **Cylinder Valves** “off” (clockwise). Turn the system **On/Off** valves “on” and “off” repeatedly until all pressure gauges read “0” psi and the Exhaust pressure gauge needles touch the silver stop.
2. Turn system **On/Off** knobs “off”.
3. Remove probes.
4. Turn **Main Meter** switch “off”. Place D.T.T.I. switch (10) in the center position.
5. Exhaust hose can be disconnected from console for storage.

## Care and Maintenance

The ACC2 should be handled with the care appropriate to any precision medical instrument. Avoid banging, dropping or mishandling the system and probes. Periodically inspect the hoses, fittings and probes for any damage or wear. Other than maintaining an adequate supply of medical grade nitrous oxide, the system is essentially maintenance free.

### ACC2 Cleaning

**NOTE:** Do not spray or pour liquids directly on the unit.

**NOTE:** The unit and/or accessories cannot be sterilized.

**WARNING:** Ensure Isopropyl Alcohol (IPA) is completely dry before operating the unit.

**CAUTION:** Avoid caustic or abrasive cleaners

#### **Guidelines**

The following guidelines are recommended for cleaning the unit. It is the user's responsibility to qualify any deviations from these processing methods.

1. Disconnect the unit or cart from the outlet before cleaning.
2. If the unit and/or accessories are contaminated with blood or other body fluids, they shall be cleaned before the contamination can dry (within two hours of contamination).
3. The outer surfaces of the unit and/or accessories shall be cleaned with 70% -90% Isopropyl alcohol (IPA) wipes for a minimum of two minutes. Do not allow fluids to enter the chassis.
4. Pay attention to all areas where fluids or soil may gather, such as under/ around the handles or any tight crevices/ grooves.
5. Dry the unit and/or accessories with a dry, white lint-free cloth.
6. Conduct a final confirmation of the cleaning process by visually inspecting the white cloth for remaining soil.
7. If soil remains on the white cloth, repeat steps 3 through 6.
8. After cleaning is complete, turn the unit on to perform Power On Self-Test (POST). If any errors are received, contact AtriCure to begin return process.

### Recommended ACC2 Probe Cleaning and Sterilization

#### EO Sterilization Parameters:

Preconditioning Set Points:

- Temperature: 55°C
- Relative humidity: 80%
- Vacuum Set Points: 1.3 psia
- Preconditioning Time: 30 minutes

Sterilization Parameters:

- Temperature: 55°C
- Relative humidity: 80%
- EO concentration: 740 mg/L
- Gas exposure time: 60 minutes

Aeration:

- Aeration time: 12 hours
- Aeration temperature: 55 ± 4°C

#### STERRAD® Sterilization Instructions:

- Place device into Tyvek® pouch and sterilize in STERRAD unit for a full STERRAD 100S, NX, or 100NX cycle using STERRAD 100S standard (short) sterilization parameters.

### Storage

Store nitrous oxide cylinders, ACC2 System, and probes at room temperature.

### Cautions

1. Nitrous oxide gas must be safely exhausted. Do not restrict or limit the gas flow in any way.
2. Always ensure that the safety latch is engaged when inserting the probe plug into the probe jack on the front of the system.
3. Use Atricure part no. F032873 exhaust tube only.
4. Do not use the system below 650 psi. Low pressure will reduce performance.
5. Do not autoclave probes or use a probe if it has been autoclaved.
6. Do not remove or install probe unless the line and exhaust pressure gauge needles touch the silver stops.

The ACC2 Probes can be ethylene oxide (EO) gas sterilized or STERRAD® sterilized.

**Initial Cleaning:**

- Rinse the device in luke warm water until all gross soils are removed.
- Wipe down probe with enzymatic detergent per detergent's manufacturer's instructions using a sponge or common wiping media.
- Rinse immediately in luke warm water (not hot) until all signs of the detergent are removed.
- Dry probe with a cloth towel or lint-free cloth.
- Personal protection equipment should be worn during these steps.

Probes must be thoroughly cleaned prior to sterilization. Always follow the instructions supplied by the manufacturer of the sterilizer. Use only valid parameters, conditions, and equipment to assure sterility.

**DO NOT AUTOCLAVE PROBES.**

**DO NOT IMMERGE PROBE AND LINES IN ANY SOLUTION.**

If probe is accidentally autoclaved or immersed—**DO NOT USE**—Return to AtriCure for servicing.

# Troubleshooting

## Troubleshooting ACC2 Cryosurgical System

### ***Probe Does Not Freeze Properly:***

**Cause**

1. On/Off valve off.
2. Cylinder valves closed.
3. Cylinder pressure below 650 psi.
4. Cylinders are cold.
5. Crimped probe or exhaust line.
6. Dirty filter in probe plug.
7. Probe not properly engaged.
8. Probe tip temperature not recording on display.
9. Leaks in serviceable area.
10. Controls not in proper position.
11. Contaminated Nitrous Oxide Gas.

**Remedy**

1. Turn valve "on".
2. Open cylinders.
3. Replace cylinders or refer to 4.
4. Warm cylinders to room temperature
5. Check line for restrictions.
6. Remove surface contaminants with nylon brush and retest probe.
7. Insure probe is properly engaged and locked in position.
8. Refer to Digital Display, Section.
9. Tighten fittings. Replace nipple on cylinder yoke or replace damaged "O" rings
10. Properly set controls.
11. Replace cylinders.

12. Console or probe malfunction.
13. Clogged filter or contamination.

### ***Probe Does Not Defrost Properly:***

#### **Cause**

1. Controls not set properly.
2. Exhaust line not connected.
3. Probe does not defrost during operation.

### ***Digital Temperature Display Malfunction:***

#### **Cause**

1. Meter not "on".
2. Open circuit in thermocouple.
3. Display out of calibration.
  
4. Batteries low

12. Connect a second probe and test.
13. Replace filter, part no. F021720.

#### **Remedy**

1. Properly set controls.
2. Ensure probe plugs are properly connected.
3. Immediately shut console off and flood tip area with sterile saline.

#### **Remedy**

1. Turn main meter switch "on".
2. Connect a second thermocoupled probe to console and test.
3. Check digital meter calibration with probe connected and system off. Submerge probe tip in ice water. Reading should be approximately 0° C. No calibration required.
4. Replace four 9 volt batteries.

## ***Control Valves Do Not Operate Properly:***

### **Cause**

1. Knob is loose on shaft.
2. Valve freezing.

### **Remedy**

1. Tighten knob with allen wrench.
2. Remove knob and tighten packing nut.  
*Do Not Over Tighten. Tighten only until leak stops.*

## ***Difficult to Insert Probe Plugs into Probe sockets or Leaks:***

### **Cause**

1. No lubricant.
2. "O" rings damaged.

### **Remedy**

1. Apply a small amount of lubricant, e.g. syringe lubricant, to the probe shank; then insert shank into the console. Twist the plug while inserting to lubricate "O" ring. *Do not place lubricant on filter or in probe opening.*
2. Replace with new "O" ring. Refer to user information for part number.

## ***Exhaust Hose Slips Off Console:***

### **Cause**

1. Exhaust line restricted.
2. External exhaust source restricted or not adjusted correctly.

### **Remedy**

1. Check exhaust line for restrictions and correct
2. Check venting for restrictions and correct adjustments.

## ***Digital Timer Display Malfunction:***

### **Cause**

1. Batteries low.
2. Controls not set properly.

### **Remedy**

1. Replace batteries (Refer to timer section).
2. Properly set controls (Refer to timer section).

## Replacement Parts for ACC2 Console, Cat. No. 3050

<u>Catalog Number</u>	<u>Description</u>
<b>F032873</b>	Exhaust Tube, 20 ft
.....	Batteries, 9 volt
Exhaust Tube, 20 ft.	"O" Ring - Upper Probe Socket
<b>F010313</b>	"O" Ring - Lower Probe Socket
.....	Probe I.D. Band
Batteries, 9 volt	22 gauge Needles – Red
<b>F020308</b>	22 gauge Needles – Yellow
.....	30 gauge Needles – Red
"O" Ring - Upper Probe Socket	30 gauge Needles – Yellow
<b>F020440</b>	Filter Cartridge – Replacement
.....	Filter "O" Ring – Replacement
"O" Ring - Lower Probe Socket	Pressure Gauge - Filter Assembly
<b>F033523</b>	30" Pressure Hose
.....	Yoke Tip Replacement
Probe I.D. Band	Caster – Braking
<b>F03039</b>	Base, Legs
.....	Mold Trim
22 gauge Needles - Red	Wrench
<b>F03041</b>	
.....	
22 gauge Needles - Yellow	
<b>F03037</b> (Optional)	
.....	
30 gauge Needles - Red	
<b>F03038</b> (Optional)	
.....	
30 gauge Needles - Yellow	
<b>F021720</b>	
.....	
Filter Cartridge - Replacement	
<b>F010924</b>	
.....	
Filter "O" Ring - Replacement	
<b>F051875</b>	
.....	
Pressure Gauge - Filter Assembly	
<b>F020370-3</b>	
.....	
30" Pressure Hose	
<b>F021837</b>	
.....	
Yoke Tip Replacement	

**F021910**

.....  
Caster - Braking  
**F052329-202**

.....  
Base, Legs  
**F030740 .**

.....  
Mold Trim  
**F021913**

.....  
Wrenc  
.....

**Note:** The only user serviceable parts contained in the ACC2 console are the timer batteries and the temperature display meter batteries.

Manufactured by:

**AtriCure Incorporated**  
7555 Innovation Way

Mason, Ohio 45040 USA  
Customer Service:  
1-866-349-2342 (toll free)  
513-755-4100

European Representative:  
**Köntges SPRL**

Avenue Hellevelt 35  
B-1180 Brussels Belgium  
Tel: +32 (0) 2 375 51 63  
FAX: +32 (0) 2 375 89 06  
e-mail: [herbert.kontges@skynet.be](mailto:herbert.kontges@skynet.be)

# AtriCure<sup>®</sup>

## frigitronics<sup>®</sup>


### Operation and Maintenance Manual


# CCS-200

Dual Probe, Non-Electric Cardiac Cryosurgical System

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