

DESCRIPTION The ATRICOR ISOlator Surgical Ablation System is composed of the Ablator and Sensing Unit (ASU-115A303-230), an ATRICOR ISOlator device, and a bipolar Transapical System. The ISOlator is an electrical instrument designed for use only with the ASU. The ISOlator is used for cardiac tissue ablation. When activated, the ASU delivers alternating RF energy to the linear electrodes on the isolated jaws of the ISOlator. The Operator controls the application of the RF energy by pressing the footswitch.




FIGURE 1 ATRICOR ISOlator ILLUSTRATION AND NOMENCLATURE There are three families of ISOlator clamps available for use with the ASU. The ATRICOR ISOlator device is configured as avascular clamp and features clamping jaws of various lengths and curvatures. There are three families of ISOlator clamps available for use with the ASU. The ATRICOR ISOlator device is configured as avascular clamp and features clamping jaws of various lengths and curvatures.

FIGURE 2 GLIDETAP™ TAPE INSTRUMENT GUIDE ILLUSTRATION AND NOMENCLATURE There are two families of GLIDETAP™ Tape Instrument Guides for use with ATRICOR ISOlator Devices. The GLIDETAP™ Tape Instrument Guide is a single patient, surgical device designed to facilitate the guidance of surgical instruments through soft tissue using general surgical techniques. The GLIDETAP™ Tape Instrument Guide (Figure 2) is designed to fit devices that use the integral attachment tip. The GLIDETAP™ Tape Instrument Guide (Figure 2) is designed to fit Instruments without integral attachment tips. NOTE: Please refer to the ATRICOR ASU and ASB3 Instructions for use for information specific to the ASU and ASB3.

POTENTIAL COMPLICATIONS DEVICES: Possible complications related to the creation of the linear lesions in cardiac tissue using a clamp-type device may include but are not limited to: Tissue Cutting, Irregular heart rhythm disturbance (atrial and/or ventricular), Postoperative embolic complication, Pericardial effusion or tamponade, Injury to the great vessels, Injury to the coronary arteries, Injury to the pulmonary veins (SA/AV Node), Injury to unintended surrounding tissue structures, including tears and punctures, Injury to myocardial perforator or septum, Injury to surrounding surrounding tissue structures, including tears and punctures, Injury to myocardial perforator or septum, Injury to surrounding surrounding tissue structures, including tears and punctures, Injury to myocardial perforator or septum, Injury to surrounding surrounding tissue structures, including tears and punctures.

FIGURE 3 ILLUSTRATION OF THE ATRICOR ISOlator INLINE SYNERGY GRIP CLAMPS. Labels include Distal Jaw, Handle, Proximal Jaw, Attachment Tip, Release Mechanism, Electrode, Connector, and ASU Alignment Arrow. Only on devices packaged with GLIDETAP™ Tape Instrument Guide.

FIGURE 4 ILLUSTRATION OF THE ATRICOR ISOlator INLINE SYNERGY GRIP CLAMPS. Labels include Lateral Jaw, Instrument Attachment Tip, Clear Ribbon, Distal Jaw, and Accessory Ribbon.

FIGURE 5 ILLUSTRATION OF THE ATRICOR ISOlator INLINE SYNERGY GRIP CLAMPS. Labels include Lateral Jaw, Instrument Attachment Tip, Clear Ribbon, Distal Jaw, and Accessory Ribbon.

FIGURE 6 ILLUSTRATION OF THE ATRICOR ISOlator INLINE SYNERGY GRIP CLAMPS. Labels include Lateral Jaw, Instrument Attachment Tip, Clear Ribbon, Distal Jaw, and Accessory Ribbon.

FIGURE 7 ILLUSTRATION OF THE ATRICOR ISOlator INLINE SYNERGY GRIP CLAMPS. Labels include Lateral Jaw, Instrument Attachment Tip, Clear Ribbon, Distal Jaw, and Accessory Ribbon.

FIGURE 8 ILLUSTRATION OF THE ATRICOR ISOlator INLINE SYNERGY GRIP CLAMPS. Labels include Lateral Jaw, Instrument Attachment Tip, Clear Ribbon, Distal Jaw, and Accessory Ribbon.

FIGURE 9 ILLUSTRATION OF THE ATRICOR ISOlator INLINE SYNERGY GRIP CLAMPS. Labels include Lateral Jaw, Instrument Attachment Tip, Clear Ribbon, Distal Jaw, and Accessory Ribbon.

FIGURE 10 ILLUSTRATION OF THE ATRICOR ISOlator INLINE SYNERGY GRIP CLAMPS. Labels include Lateral Jaw, Instrument Attachment Tip, Clear Ribbon, Distal Jaw, and Accessory Ribbon.

FIGURE 11 ILLUSTRATION OF THE ATRICOR ISOlator INLINE SYNERGY GRIP CLAMPS. Labels include Lateral Jaw, Instrument Attachment Tip, Clear Ribbon, Distal Jaw, and Accessory Ribbon.

FIGURE 12 ILLUSTRATION OF THE ATRICOR ISOlator INLINE SYNERGY GRIP CLAMPS. Labels include Lateral Jaw, Instrument Attachment Tip, Clear Ribbon, Distal Jaw, and Accessory Ribbon.

FIGURE 13 ILLUSTRATION OF THE ATRICOR ISOlator INLINE SYNERGY GRIP CLAMPS. Labels include Lateral Jaw, Instrument Attachment Tip, Clear Ribbon, Distal Jaw, and Accessory Ribbon.

FIGURE 14 ILLUSTRATION OF THE ATRICOR ISOlator INLINE SYNERGY GRIP CLAMPS. Labels include Lateral Jaw, Instrument Attachment Tip, Clear Ribbon, Distal Jaw, and Accessory Ribbon.

FIGURE 15 ILLUSTRATION OF THE ATRICOR ISOlator INLINE SYNERGY GRIP CLAMPS. Labels include Lateral Jaw, Instrument Attachment Tip, Clear Ribbon, Distal Jaw, and Accessory Ribbon.

FIGURE 16 ILLUSTRATION OF THE ATRICOR ISOlator INLINE SYNERGY GRIP CLAMPS. Labels include Lateral Jaw, Instrument Attachment Tip, Clear Ribbon, Distal Jaw, and Accessory Ribbon.





