

INSTRUCTIONS FOR COMPLETION:

The following information should be completed by the healthcare institution/provider on the front of the implant card (See **FIGURE 1**).

	Patient Identification
	Procedure Date
	Name and Address of the Implanting Healthcare Institution/Provider

Once the implant card is populated with the information above, peel the card from the backer, fold along the perforation, and bond it together. Give both the implant card leaflet and implant card to the patient.

PATIENT INFORMATION:

The implant card serves as a record of the Atrial Closure implant you have received. The information on the implant card helps communicate information about this implant. It can be presented to physicians or airport security. The implant card should be presented anytime you are subject to a Magnetic Resonance Imaging (MRI) screening, as the information helps ensure the MRI does not interfere with the implanted device.

In the event you need to get a replacement card for your implant, you will need to have the UDI and lot information for your device off the rear of the card (See **FIGURE 2**). Contact the manufacturer to request a replacement card.

In the unlikely event of any issues with the implant, the information on the implant card will ensure that your physician can inform you if you are impacted.

Device lifetime is equal to patient lifetime.

MR CONDITIONAL

Non-clinical testing demonstrated that the LAA Exclusion System Clip is MR Conditional. A patient with this device can be scanned safely in an MR system immediately after placement under the following conditions:

- Static magnetic field of 1.5-Tesla and 3-Tesla, only.
- Maximum spatial gradient magnetic field of 4,000-gauss/cm (40-T/m)(extrapolated) or less.
- Maximum MR system reported, whole body averaged specific absorption rate (SAR) of 4-W/kg for 15 minutes of scanning (i.e., per pulse sequence) in the First Level Controlled Operating Mode of operation for the MR system.
- The scan conditions defined for the LAA Exclusion System Clip are expected to produce a maximum temperature rise of 2.9°C (5.22°F) for Gillinov-Cosgrove Clip (ACH1, ACH2, PRO1, PRO2 devices) and 3.1°C (5.58°F) for V Clip (FLEX-V, PRO-V devices) after 15-minutes of continuous scanning (i.e., per pulse sequence).

DEVICE CONTAINS		SYMBOLS GLOSSARY			
CGC CLIP	V CLIP		Device Name		Unique Device Identifier
Titanium Grade 2	Titanium Grade 5		Batch Code		Information Website for Patients
Polyurethane	Polyethylene Terephthalate		Manufacturer		MR Conditional
Nitinol	Titanium Dioxide		Meets the Medical Device Regulation (MDR) 2017/745		Authorized representative in the European Community
Polyethylene Terephthalate					
Titanium Dioxide					

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