# AtriCure®

# CE

# 2797

IMPLANT CARD LEAFLET

**MD**

**FIGURE 1**

	_____
	_____
	_____
	_____
	_____
	https://www.atricure.com/international/instructions-for-use
<b>AtriCure®</b>	PL-001.H

**FIGURE 2**



<b>en</b> Atrial closure / <b>bg</b> Затваряне на предсърдие / <b>cs</b> Uzávěr srdeční síně / <b>da</b> Lukning af forkamre / <b>de</b> Vorhof-Verschluss / <b>el</b> Κοιλιακή σύγκλιση / <b>es</b> Cierre de la aurícula / <b>et</b> Kodade sulgemisseade / <b>fi</b> Eteläisen sulkuale / <b>fr</b> Fermeture auriculaire / <b>hr</b> Zatvaranje atrija / <b>hu</b> Pitvari lezárás / <b>is</b> Gátta lokun / <b>it</b> Chiusura atriale / <b>lt</b> Prieširdžio uždarymas / <b>lv</b> Priekškambara slēgums / <b>nl</b> Atriale sluiting / <b>no</b> Lukkemekanisme til atrium / <b>pl</b> Zamknięcie przedsionka / <b>pt</b> Encerramento auricular / <b>ro</b> Închidere atrială / <b>ru</b> Закрытие предсердия / <b>sk</b> Uzatváranie predsienie / <b>sl</b> Atrijsko zapiralo / <b>sr</b> Zatvaranje pretkomore srca / <b>sv</b> Stängning av hjärtformak / <b>tr</b> Atrial Kapatma	
<b>MD</b> <b>ACH240</b> <b>AtriClip ACH240 40mm</b>	
<b>LOT</b> <b>1234567890</b>	<b>UDI</b>
 AtriCure, Inc. 7555 Innovation Way Mason, Ohio 45040 USA +1 866 349 2342 +1 513 755 4100	UDI-DI:(01) <b>10840143919695</b> 

IMPLANT CARD LEAFLET

**en**

### INSTRUCTIONS FOR COMPLETION:

The following information should be completed by the healthcare institution/provider on the front of the implant card (**See FIGURE 1**).

	Patient Identification
	Procedure Date
	Name and Address of the Implanting Healthcare Institution/Provider

Once the implant card is populated with the information above, peel the card from the backer, fold along the perforation, and bond it together. Give both the implant card leaflet and implant card to the patient.

### PATIENT INFORMATION:

The implant card serves as a record of the Atrial Closure implant you have received. The information on the implant card helps communicate information about this implant. It can be presented to physicians or airport security. The implant card should be presented anytime you are subject to a Magnetic Resonance Imaging (MRI) screening, as the information helps ensure the MRI does not interfere with the implanted device.

In the event you need to get a replacement card for your implant, you will need to have the UDI and lot information for your device off the rear of the card (**See FIGURE 2**). Contact the manufacturer to request a replacement card.

In the unlikely event of any issues with the implant, the information on the implant card will ensure that your physician can inform you if you are impacted.

Device lifetime is equal to patient lifetime.

### MR CONDITIONAL

The LAA Exclusion System Clip is MR Conditional. A patient with the LAA Exclusion System Clip may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.

- Static Magnetic Field Strength (T): 1.5-Tesla or 3.0-Tesla
- Static Magnetic Field Orientation: Horizontal

- Maximum Spatial Field Gradient (T/m and gauss/cm): 40-T/m (4,000-gauss/cm)
- RF Excitation Polarization: Circularly Polarized (CP) (i.e., Quadrature-Transmission)
- Transmit RF Coil: Any transmit RF coil may be used
- Receive RF Coil: Any receive-only RF coil may be used
- MR System Operating Mode: Normal Operating Mode
- Maximum Whole Body Averaged SAR (W/kg): 2-W/kg (Normal Operating Mode)
- Scan Duration and Wait Time: Whole body averaged SAR of 2-W/kg for 60 minutes of continuous RF exposure (i.e., per pulse sequence or back-to-back sequences/series without breaks)
- MR Image Artifact: The presence of the implant may produce an MR image artifact. Imaging protocol modifications may be necessary to compensate for the MR image artifact.

		
ATRICLIP: ACH2, PRO1, PRO2	V-Clip: PROV, ACHV	ATRICLIP MINI: ACHM, PROM

DEVICE CONTAINS			SYMBOLS GLOSSARY	
<b>ATRICLIP</b>	<b>V-CLIP</b>	<b>ATRICLIP MINI</b>	<b>MD</b>	<b>UDI</b>
Titanium Grade 2	Titanium Grade 5	Titanium Ti64	Device Name	Unique Device Identifier
Polyurthane	Polyethylene Terephthalate	Polyethylene Terephthalate	<b>LOT</b>	 Information Website for Patients
Nitinol	Titanium Dioxide	Titanium Dioxide	 Manufacturer	 MR Conditional
Polyethylene Terephthalate		Nitinol	 Meets the Medical Device Regulation (MDR) 2017/745	 Authorized representative in the European Community
Titanium Dioxide				

Not all products/models may be available in your specific country. Contact your local AtriCure representative to check availability.

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AtriCure Inc.
AtriCure Europe B.V.