
# LPK-302

## External Graphics Display (EGD) Software

*for use with the CSK-310 RF Generator*

## OPERATOR'S MANUAL

AtriCure, Inc.  
7555 Innovation Way  
Mason, Ohio 45040, USA  
1.866.349.2342  
www.atricure.com

AtriCure Europe B.V.  
De entree 260  
1101 EE Amsterdam  
NL  
+31 20-7005560  
ear@atricure.com

**TABLE OF CONTENTS**

**Description of Contents**.....2  
**Required But Not Provided** .....2  
**Indications for Use**.....2  
**Intended Purpose** .....2  
**Intended User** .....2  
**Target Patient Population** .....2  
**Clinical Benefit**.....2  
**Serious Incident Statement**.....2  
**Contraindications** .....2  
**Data Archival** .....2  
**Warnings And Cautions** .....2  
**Symbol Glossary** .....2  
**Instructions for Use** .....3  
**1. Device Description / System Overview**.....3  
    1.1. RF EGD Software Installation..... 3  
    1.2. Software Interface/Features and Functions..... 3  
        1.2.1. RF\_Display software ..... 3  
        1.2.2. Selecting Procedure ..... 4  
        1.2.3. Selecting Procedure ..... 4  
        1.2.4. Selecting Tools ..... 5  
        1.2.5. Selecting Windows..... 5  
        1.2.6. Selecting Language..... 5  
        1.2.7. Closing the Software ..... 6  
    1.3. Connecting RF Generator to the Laptop..... 6  
    1.4. Operating Modes..... 6  
**2. Technical Specifications**.....6  
    2.1. RF EGD Software Version / Device Classification ..... 6  
    2.2. Laptop Specifications ..... 6  
    2.3. Expected Lifetime ..... 6  
    2.4. Environmental Specifications..... 6  
**3. System Disassembly After Use** .....6  
**4. Cleaning and Preventive Maintenance** .....6  
    4.1. Cleaning and Disinfecting Instructions..... 6  
    4.2. Product Support..... 6  
    4.3. Disposal ..... 6  
**SUMMARY OF SAFETY AND CLINICAL PERFORMANCE**     **6**  
**5. Troubleshooting** .....7  
**Disclaimer** .....7

**INDICATIONS FOR USE**

The External Graphics Display (EGD) software is an optional accessory indicated to be used in conjunction with the CSK-310 Generator System to display the energy delivered during each ablation.

**INTENDED PURPOSE**

The External Graphics Display (EGD) software is an optional accessory intended to be used in conjunction with the CSK-310 Generator System to display the energy delivered during each ablation.

**INTENDED USER**

Licensed medical doctors who perform cardiac and/or thoracic procedures using AtriCure ablation handpieces.

**TARGET PATIENT POPULATION**

Adult patients with arrhythmias including atrial fibrillation.

**CLINICAL BENEFIT**

To achieve the clinical benefits of the attached AtriCure ablation handpieces.

**SERIOUS INCIDENT STATEMENT**

Any serious incident that has occurred in relation to this device should be reported to AtriCure and the competent authority of the Member State in which the user and/or patient is located.

**CONTRAINDICATIONS**

None

**DATA ARCHIVAL**

**Important:** The data displayed on the Generator/software screen can be exported as a Comma Separated Value (CSV) file for archival purposes.

**WARNINGS AND CAUTIONS**

 <b>WARNINGS</b> 
<p>Read all instructions carefully prior to using the device. Failure to do so may result in injury and/or improper device function.</p> <p>Information displayed on the LPK-302 software should not be used alone to make treatment decisions, making such decisions based on the display alone could lead to patient injury.</p>

 **CAUTIONS**

- The software should be used by physicians or users trained in the proper use of Operating Room (OR) equipment.
- Use only with the CS-3000 RF Generator do not use with any other systems, to prevent injury and/or equipment damage.
- Wipe the non-sterile USB Stick per the Cleaning Instructions mentioned in Section 4.1 before installing the software on the laptop detailed in Section 1.1. Retain the USB outside the OR.
- The software installed on a laptop serves only as a secondary display projecting information from the Generator in a different format. The controls within the software/laptop do not operate the Generator.
- No modification of this software/equipment is allowed.
- Do not spray or pour liquids directly on or allow fluids to enter the USB Stick containing the Software Installer/laptop.

**SYMBOL GLOSSARY**

	Conforms to the requirements of the European Directives and Regulations		Waste Electrical and Electronic Equipment (WEEE)
	Lot Number		Medical Device
	Authorized Representative in the European Community		Importer
	Caution		Consult Instructions for Use



**INSTRUCTIONS FOR USE** en

**DESCRIPTION OF CONTENTS**

- One (1) USB Stick containing Software Installer (provided non-sterile)
- One (1) Operator's Manual

**REQUIRED BUT NOT PROVIDED**

- Laptop (with minimum requirements per Section 2.2) and its accessories
- nContact CSK-310 RF Generator (also known as CS-3000 when referencing generator controls) and components
- Data cable (standard, commercially available, USB 2.0 A/B Cable with Ferrite Chokes (Male/Male))

	Manufacturer		Date of Manufacture
	Disconnected (Data cable connection between the Generator and the laptop)		Navigation buttons to review previous/next ablation(s)
	Shortcut to RF_Display application (created and located on laptop desktop)		Generator Connections Status Indicator (Green or Red)
	Model Number		Non-Sterile
W	Watts	S	Seconds
$\Omega$	Ohms		Keep dry
	Unique Device Identifier		Catalog Number
	Does Not Contain Latex		Does Not Contain Phthalates
 <p>Humidity Limits</p>		 <p>Transit Temperature Limits</p>	

## INSTRUCTIONS FOR USE

### 1. DEVICE DESCRIPTION / SYSTEM OVERVIEW


The RF EGD software installed on a laptop is an optional accessory that provides ablation information such as a graphical display of Power and Impedance and total energy delivered (Max Joules).

#### 1.1. RF EGD Software Installation

**NOTE:** The USB has been Quality-checked by the manufacturer for virus/malware. Do not use if seal is broken.

1. Remove the non-sterile USB Stick from the shipper box and check for any signs of material degradation, such as corrosion, bending of the connection area, or any cracks. If any signs of degradation, do not attempt software installation and contact AtriCure Customer Service (refer to Section 4.2). Wipe the non-sterile USB Stick per the Cleaning Instructions mentioned in Section 4.1 before proceeding with the software installation.
2. Plug in the USB Stick into any one of the laptop's USB ports and open the contents as shown below. Double-click the executable file in Figure 1.

**FIGURE 1: PASSWORD-PROTECTED FOLDER WITHIN THE USB STICK CONTAINING THE SOFTWARE INSTALLER**

3. Call AtriCure Customer Service (refer to Section 4.2) to obtain the password. Enter the password when prompted and click Login as shown in Figure 2.

**FIGURE 2: ENTERING THE PASSWORD TO OPEN THE SOFTWARE INSTALLER**


4. Right-click the RF\_Display folder and select Copy. Paste the RF\_Display folder on the desktop and double-click the folder to open it. Double-click the RF\_Display\_Setup executable file and follow the instructions for installation.
5. Check the type of data port (USB) on the backside of the RF Generator:
  - a) For the USB data port as shown in Figure 12, double-click Attachment 2, and run the executable present within. Click Extract and follow instructions to install.
6. Click "Log-off" in the top left corner and then close the SanDiskSecureAccess window. Right-click the USB icon in the bottom right corner of the laptop and select Eject. Remove the USB Stick from the laptop when the Safe to remove the hardware prompt appears on the screen. Retain the USB outside the Operating Room (OR).
7. Restart the laptop before using the application for the first time.

**NOTE:** One-time software installation is necessary for every laptop intended for use with the CS-3000 RF Generator. It is essential to restart the laptop before using the software for the first time. Retain the USB Stick and this Operator Manual for accessing the password-protected folder within the USB Stick to install this software on additional laptops as required.

**NOTE:** The software maximizes and occupies the complete screen of the laptop. For the best visual experience, the recommended screen resolution is 1280 x 720 pixels.

### 1.2. Software Interface/Features and Functions

**FIGURE 3: SOFTWARE INTERFACE FEATURES/MENU**


#### 1.2.1. RF\_Display software has five different features/menus in the software (refer Figure 3):


1. Procedure:
  - a) New: Opens a new window/session for capturing data from the Generator and plotting the Power-Impedance over Time graph and energy delivered in the bar graph.
  - b) Open: Opens a window for reviewing the previously performed ablations.
2. Tools:
  - a) Configure Port: User settings that help establish communication between the Generator and the laptop.
3. Windows:
  - a) New Window: Opens a new window/session for capturing data from the Generator and plotting the Power-Impedance over Time graph and energy delivered in the bar graph.
  - b) Cascade: Arranges all open windows in a cascading format.
  - c) Tile Vertical: Arranges all open windows vertically.
  - d) Tile Horizontal: Arranges all open windows horizontally.
  - e) Close All: Closes all open windows but does not shut down the software.
4. Language:
  - a) User can change all the on-screen text (within the software) to any language by selecting country icon.
5. Help:
  - a) Displays Application name, Software version, Manufacturer name and address, Contact details, and Device Identifier for the software.

#### 1.2.2. Selecting Procedure > New opens the following window (refer Figure 4 and Figure 5):

**FIGURE 4: OPTION FOR OPENING A NEW PROCEDURE WINDOW**

**FIGURE 5: SOFTWARE INTERFACE LAYOUT**

1. Device Length Indicator	8. Impedance Graph (yellow)
2. Lesion Time Indicator	9. Power Graph (blue)
3. Lesion Power Indicator	10. Disconnected (red)
4. Reason for End of Lesion	11. Review Lesion(s) (< >)
5. Lesion Number Indicator (User input)	12. Reset Button (resets graph)
6. Communication Status Indicator	13. Energy Delivered Bar Graph
7. Lesion Number Indicator	

#### Device, Time, and Set Power Settings

Device settings (power/time) are pre-programmed in the EPI-Sense® coagulation device and are recognized by the Generator upon plugging in the device. These default settings appear on the software (under Time and Set Power, respectively). Device in the software displays the device length.

**NOTE:** Verify the Power and Time values are displayed in the Graphical window before starting an ablation.

#### Reason for End of Lesion

The Generator continues to deliver RF energy for ablation unless one of the three criteria is satisfied. These reasons will be displayed in the “Reason for End of Lesion” window [#4 in Figure 5]:

1. Time Limit Reached: Time limit specified under Time in the software has been reached.
2. Impedance limits reached: If tissue impedance falls below 30 Ω or increases beyond 500 Ω.
3. User terminated: Occurs when the user presses the RF ON/OFF button or foot pedal at any time during an ongoing ablation.

#### Communication Status Indicator

The Communication Status indicator provides an easy visual aid regarding the communication status between the Generator and the laptop. This indicator continues to blink red if the communication is not established (refer to Section 5: Troubleshooting for potential causes and actions for establishing communication between the Generator and the laptop). The indicator turns to solid green as soon as the communication gets established.

#### Lesion

It is a numeric field capable of accepting up to six digits as user input. Users can enter and change the lesion number for each lesion. This number also appears in the bottom left corner of the software’s graph window as soon as the ablation ends and during the review of ablations performed.

#### Maximum Percent Joules of Energy Delivered Bar Graph

The calculation for the energy delivered by the Generator displayed as a bar graph (as percent Maximum Joules) is as follows:

$$\% \text{ Max Joules} = \frac{(\text{Achieved Power} \times \text{Achieved Time}) \times 100}{(\text{Set Power} \times \text{Set Time})}$$

**Note:** Power output of the generator must be within +/-20% of the setpoint.

#### Power-Impedance Over Time Graph

The power delivered by the Generator and the tissue impedance, both trending over time, are displayed in the software in a graph with blue and yellow colors, respectively. The plotting stops when the data cable between the Generator and the laptop gets disconnected until the connection reestablishes.

The **Disconnected** icon  will blink for the entire ablation from the time of disconnection. Even though the real-time plotting may stop within the Software, the Generator continues to provide RF energy to the tissue as intended.

The user may review the ablation data by selecting Windows > Close All to close the current ablation window and select Procedure > Open > Procedure file > OK or use Browse to open the procedure file.

The data between the time points during which the Generator disconnected from the laptop displays in red. During such time points, Power and Impedance data will not get saved on the laptop, and values for both Power and Impedance are shown as zero even though the Generator continues to provide RF energy as intended.


#### Reset button

The user may press the Reset button to clear both the plotted Power-Impedance over Time graph and Percent Maximum Joules bar graph before beginning the next ablation. Both these graphs automatically clear when the new ablation begins.


#### 1.2.3. Selecting Procedure > Open (refer Figure 6 and Figure 7):

Previous ablations can be reviewed either by pressing <> in the current Procedure window or by using Procedure > Open > Procedure file (Procedure date, time, and language of on-screen text).

**FIGURE 6: OPTION FOR OPENING SPECIFIC ABLATION FOR REVIEW**


**FIGURE 7: SELECTING PREVIOUSLY CONDUCTED ABLATION FOR REVIEW (IMAGE FOR REFERENCE ONLY)**


1.2.4. Selecting Tools > Configure Port opens the following window (refer Figure 8 and Figure 9):

**FIGURE 8: OPTION FOR CONFIGURING COMMUNICATION PORT SETTINGS IN THE SOFTWARE**

**FIGURE 9: OPTIONS FOR CONFIGURING COMMUNICATION PORT SETTINGS IN THE SOFTWARE**


If the Communication Status indicator is blinking red in the Procedure window, the user should perform any of the following actions:

- a) Ensuring the data cable is connected to both the Generator and the laptop
- b) Connecting the data cable to a different USB port on the laptop
- c) Selecting the other option under the Stop Bits dropdown settings.

Refer to Section 5: Troubleshooting for potential causes and actions for establishing communication between the Generator and the laptop).

1.2.5. Selecting Windows opens the following window (refer Figure 10):

**FIGURE 10: OPTIONS FOR ARRANGING WINDOWS IN THE SOFTWARE**

Users may arrange various open windows within the software in either a cascading, horizontal, or vertical format. Any of these windows may be closed individually or together using the Close All option.

1.2.6. Selecting Language (refer Figure 11):

**FIGURE 11: OPTIONS FOR SELECTING LANGUAGE IN THE SOFTWARE FOR ALL ON-SCREEN TEXT TRANSLATION**

Users may change the language of all the on-screen text to any of the languages provided. The default language of the software is English.

The software automatically closes and restarts with all the on-screen text translated into the selected language upon selecting a different language other than the current language.

If the user decides to cancel the translation, the software does not close and goes back to the home screen (as shown under Section 1.2). A prompt alerts the user if the user selects the same language for translation. Flags of countries are available alongside language names for additional user convenience.

**NOTE:** Language selection is retained even upon closure of the software.

### 1.2.7. Closing the Software

The user must close the software by pressing X in the top right corner of the application at the end of each procedure.

**NOTE:** Application cannot be closed during an ablation. An alert requests the user to complete the current (or manually terminate) the ablation. Completion of an ablation is essential to allow the saving of data for future reference.

### 1.3. Connecting RF Generator to the Laptop

Data cable connections on the rear side of the Generator are per the illustration shown in Figure 12.

**FIGURE 12: nContact CS-3000 RF GENERATOR REAR PANEL WITH DATA PORT (USB) (IMAGE FOR REFERENCE ONLY)**

1. USB Data Port

**Caution:** Please see CS-3000 RF Generator Unit IFU for additional system information.

**NOTE:** Set up the CS-3000 RF Generator with all the cables connected (power cord and USB-USB data cable between the Generator and the laptop) before starting the procedure.

### 1.4. Operating Modes

The software works only when the Generator is in the READY POWER CONTROL mode and RF ON POWER CONTROL modes.

## 2. TECHNICAL SPECIFICATIONS

### 2.1. RF EGD Software Version / Device Classification

- Device Classification: Class IIa
- Software Installer version: 3.3.1

### 2.2. Laptop Specifications

Operating System	Windows 10 (Home, Enterprise, and Pro versions only) or newer
RAM	4 GB or greater
Hard Disk Drive	32 GB hard drive or 32 GB SSD or greater
Processor Speed	At least 1.5 GHz
Screen Resolution	1280 x 720 resolution monitor
Screen Size	At least 11 inches
Size	14 in (356 mm) - (L) x 10 in (254 mm) - (W) x 1.00 in (26 mm) - (H) maximum
USB Ports	At least 2 (USB 2.0 and 3.0 enabled/compatible)
.NET Framework	4.5.1 minimum

### 2.3. Expected Lifetime

AtriCure has defined the Expected Lifetime of LPK-302 to be 5 years.

### 2.4. Environmental Specifications

Operating Conditions	
Temperature	0°C to 35°C, 32°F to 95°F
Humidity	15% RH to 85% RH, non-condensing
Storage Conditions	
Temperature	-10°C to 70°C, 14°F to 158°F
Humidity	15% RH to 85% RH, non-condensing
Transit Conditions	
Temperature	-10°C to 70°C, 14°F to 158°F
Humidity	15% RH to 85% RH, non-condensing

## 3. SYSTEM DISASSEMBLY AFTER USE

- After completing all ablations for that procedure, close the software by pressing the X in the top right corner of the application.
- Unplug the data cable from the back panel of the Generator.
- Turn off the laptop and unplug the power cord.

## 4. CLEANING AND PREVENTIVE MAINTENANCE

### 4.1. Cleaning and Disinfecting Instructions

- If the USB Stick is contaminated with blood or other body fluids, clean them before the contamination can dry.
- Dry the USB Stick with a dry, soft, clean, white lint-free cloth.
- Check for signs of material degradation, such as corrosion, bending of the connection area, or any cracks. If any signs of degradation, discard the USB Stick without attempting any software installation and contact AtriCure Customer Service (refer to Section 4.2).

### 4.2. Product Support

No preventive maintenance for the USB Stick containing the Software Installer is required. AtriCure is not responsible for any preventive maintenance for the user's laptop.

For servicing issues and to get password information, contact AtriCure, Inc. at:

**European Union (EU)**  
**AtriCure Europe BV**  
 SalesSupportEU@AtriCure.com  
 Customer Service/ Product Inquiries  
 Telephone: +31 20 7005560

### 4.3. Disposal

The USB Stick does not contain hazardous substances. Follow local governing ordinances and recycling plans regarding disposal or recycling of device components.

## SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

A summary of the safety and clinical performance of the device can be found in the European database on medical devices (EUDAMED) at <https://ec.europa.eu/tools/eudamed> by entering the Basic UDI-DI associated with the device.

Product Code(s)	Basic UDI-DI
LPK-302	0840143900000000000009ZT

## 5. TROUBLESHOOTING

Problem	Possible Cause	Action
Software installation is not successful	<ul style="list-style-type: none"> <li>The laptop configuration may not meet the recommended minimum specifications.</li> </ul>	<ul style="list-style-type: none"> <li>Verify laptop compatibility by referring to Section 2.2 for Laptop Specifications.</li> <li>Call AtriCure service if problem is not resolved.</li> </ul>
Communication indicator is flashing red	<ul style="list-style-type: none"> <li>The data cable may not have been connected between the Generator and the laptop.</li> <li>The data cable may have been disconnected either from the Generator or the laptop.</li> </ul>	<ul style="list-style-type: none"> <li>Ensure the data cable is connected securely between the Generator and the laptop.</li> </ul>
	<ul style="list-style-type: none"> <li>The data cable may be damaged.</li> </ul>	<ul style="list-style-type: none"> <li>Replace the data cable, connect the cable securely, and check if the indicator turns solid green.</li> </ul>
	<ul style="list-style-type: none"> <li>The laptop USB port may not be functioning correctly.</li> </ul>	<ul style="list-style-type: none"> <li>Remove the data cable from the laptop. Turn OFF and turn ON the Generator, reconnect the data cable to a different laptop USB port, and restart the RF_Display Software.</li> <li>If the problem continues, click Tools &gt; Configure Port, select the other option for Stop Bits, and click OK.</li> <li>Shut down and restart the RF_Display Software.</li> <li>If the problem continues, shut down the software, disconnect the data cable, and use the RF Generator display.</li> </ul>
	<ul style="list-style-type: none"> <li>The Generator may be in STANDBY mode.</li> </ul>	<ul style="list-style-type: none"> <li>Press the MODE button present on the Generator to transition from STANDBY to READY POWER CONTROL mode.</li> </ul>
Preset device settings displayed incorrectly	<ul style="list-style-type: none"> <li>Data packet from RF generator may be incomplete or corrupted</li> </ul>	<ul style="list-style-type: none"> <li>If 0W and 0s are displayed under Power and Time before beginning an ablation, the application should be closed using X (top right corner of the application), and restart the RF_Display software.</li> <li>If the problem continues, application should be closed and the generator should be power cycled. Press the MODE button to transition generator to the READY POWER CONTROL mode and restart the RF_Display software.</li> <li>If the problem continues, shut down the software, disconnect the data cable, and use the RF Generator display.</li> </ul>

“locomp Evaluation” message box pops when Procedure>New is selected	<ul style="list-style-type: none"> <li>Redundant and outdated locomp Evaluation Copy message</li> </ul>	<ul style="list-style-type: none"> <li>Move the green message box to a corner of the laptop monitor, without clicking OK or closing it with X in the top right corner.</li> </ul>
The graph does not populate during ablation	<ul style="list-style-type: none"> <li>Communication between the Generator and the laptop may be broken.</li> </ul>	<ul style="list-style-type: none"> <li>Ensure the communication indicator is solid green (if it is flashing red, follow the steps above). If the indicator is green and the grid appears to shake, continue the current ablation until the end, close the window using Windows &gt; Close All and then select Procedure &gt; New before beginning the next ablation.</li> <li>If communication was lost during ablation,  will continue to blink. Ablation data may be reviewed using Select Windows &gt; Close All, and then select Procedure &gt; Open &gt; Procedure file &gt; OK or use Browse to open the procedure file.</li> <li>If the problem continues, shut down and restart the RF_Display Software.</li> <li>If the problem continues, shut down the software, disconnect the data cable, and use the RF Generator display.</li> </ul>
Ablation data for the current case cannot be accessed	<ul style="list-style-type: none"> <li>The data cable may have been disconnected either from the Generator or the laptop during ablation.</li> </ul>	<ul style="list-style-type: none"> <li>Ensure the data cable is connected securely between the Generator and the laptop throughout ablation.</li> <li>The file path for the data log may be different than the default path.</li> </ul>
	<ul style="list-style-type: none"> <li>Attempt to close the software may have been made during ablation.</li> </ul>	<ul style="list-style-type: none"> <li>Ensure no attempt is made to close the RF_Display Software during ablation. The application should be closed using X (top right corner of the application) only after completing (or manually terminating) the current ablation to prevent loss of data.</li> </ul>
Device not recognized	<ul style="list-style-type: none"> <li>Broken/bent pins within the device connector.</li> <li>The coagulation device identification code may be out of acceptable range</li> </ul>	<ul style="list-style-type: none"> <li>If <b>Invalid device</b> or <b>None</b> is displayed under Device on the laptop screen, confirm the same is displayed on the Generator screen. If yes, follow troubleshooting information in the CS-3000 Operator Manual.</li> <li>If the problem continues, replace the coagulation device with another sterile AtriCure coagulation device.</li> <li>If the problem continues, shut down the software, disconnect the data cable, and use the RF Generator display.</li> </ul>

**DISCLAIMER**

Users assume responsibility for approving the acceptable condition of this product before it is used, and for ensuring that the product is only used in the manner described in these Instructions for Use.

Under no circumstances will AtriCure, Inc. be responsible for any incidental, special, or consequential loss, damage, or expense, which is the result of the deliberate misuse of this product, including any loss, damage, or expense which is related to personal injury or damage to property.

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För serviceproblem och för att få lösenordsinformation, kontakta AtriCure, Inc. på: