AtriClip PRO-Mini® LAA Exclusion System INSTRUCTIONS FOR USE

PROM35, PROM40, PROM45, PROM50
MD

FIGURE 1


FIGURE 2

FIGURE 3



FIGURE 4

FIGURE 5




FIGURE 6

FIGURE 7

INSTRUCTIONS FOR USE

en

AtriClip® LAA Exclusion System

INTENDED PURPOSE

The AtriClip LAA Exclusion System facilitates delivery and placement of AtriClip device for exclusion of the heart's left atrial appendage.

INDICATION FOR USE

The AtriClip LAA Exclusion System is indicated for use in patients at high risk of thromboembolism for whom left atrial appendage exclusion is warranted.

CONTRAINDICATIONS

- Do not use this device as a contraceptive tubal occlusion device.
- Do not use this device if the patient has a known allergy to Nitinol (nickel titanium alloy).
- Do not use this device if evidence of systemic infection, bacterial endocarditis, or in presence of infected operating field.

INTENDED USER

Licensed medical doctors who perform cardiac and/or thoracic procedures.

INTENDED PATIENT POPULATION

Patients at high risk of thromboembolism who are anatomically eligible for left atrial appendage exclusion.

INTENDED CLINICAL BENEFIT

Elimination of left atrial appendage, a source of thrombus, resulting in reduction in thromboembolic events.

SYSTEM DESCRIPTION

The AtriClip PRO-Mini LAA Exclusion System contains the AtriClip Mini (Clip) for exclusion of the heart's left atrial appendage (LAA).

The AtriClip PRO-Mini LAA Exclusion System is a delivery and deployment device preloaded with an AtriClip Mini. The AtriClip Mini is preloaded on a disposable Clip applicator. The AtriClip Mini is a permanent implant; device lifetime is equal to patient lifetime. The Clip was determined to be "MR Conditional" per the requirements of standard ASTM F2503-23e1.

The AtriClip PRO-Mini LAA Exclusion System is used to deliver a preloaded Clip to the target LAA site. The Clip is a sterile, permanent implant composed of Titanium (Ti64) beams, nitinol springs and covered in a knit-braided Polyethylene Terephthalate fabric that contains a small fraction of titanium dioxide. The AtriClip PRO-Mini LAA Exclusion System is not made with natural rubber latex and does not contain phthalates. Detailed materials information for implanted AtriClip Mini sizes 35 mm to 50 mm is below:

Material	Mass (g)	CAS #
Titanium Ti64	0.76-1.07	Titanium, 7440-32-6 Aluminum, 7429-90-5 Vanadium, 7440-62-2
Polyethylene Terephthalate	0.14-0.18	25038-59-9
Nitinol	0.17-0.18	Nickel, 7440-02-0 Titanium, 7440-32-6
Titanium Dioxide	0.001	13463-67-7

ENVIRONMENTAL SPECIFICATIONS

Operating	Storage	Transit
Temperature: > 20°C (68°F)	Temperature: -29°C/ -20°F to 60°C/ 140°F	Temperature: -29°C/ -20°F to 60°C/ 140°F
Relative Humidity: N/A	Relative Humidity: 15% to 85%	Relative Humidity: 30% to 85%
Atmospheric Pressure: N/A	Atmospheric Pressure: N/A	Atmospheric Pressure: N/A

PACKAGE CONTENTS

- One (1) AtriClip PRO-Mini LAA Exclusion System
- One (1) Implant Card and (1) Implant Card Leaflet

System is supplied STERILE and NON-PYROGENIC in an unopened, undamaged package. For single use only. Do not re-sterilize. Do not re-use.

SYSTEM ACCESSORIES

Other devices, not included with the System, may be used in conjunction with the AtriClip PRO-Mini LAA Exclusion System. These may include but are not limited to the following:

- Selection Guide (CGG100) (Guide)—Packaged Separately
- Minimum 12 mm port

ATRICLIP PRO-MINI LAA EXCLUSION SYSTEM

NOMENCLATURE (SEE FIGURE 1)

- | | |
|-----------------------------|---------------------------|
| [1] AtriClip Mini (Clip) | [7] Articulation Lock |
| [2] Clip Opening Jaws | [8] Lever Release Trigger |
| [3] End Effector | [9] Activation Lever |
| [4] Shaft | [10] Handle |
| [5] Left/Right Articulation | [11] Deployment Tab |
| [6] Up/Down Articulation | |

⚠️ WARNINGS ⚠️

Read all instructions carefully for the AtriClip LAA Exclusion System before use and use the device only as intended. Use of the AtriClip LAA Exclusion System should be limited to properly trained and qualified medical personnel. Improper use of this system may lead to device malfunction, failure to provide intended therapy, and/or serious injury to user or patient.

The safety and effectiveness of this device in atrial rhythm control management, either alone or in combination with ablative treatment, has not been established.

DO NOT RESTERILIZE. The AtriClip LAA Exclusion System is provided STERILE and is intended for SINGLE use only. Re-sterilization may cause loss of function or injury to patient.

Evaluate if thrombus is present in LAA. Management of thrombus is dependent on surgeon's standard of care. It is not recommended to place Clip on LAA if there is evidence of thrombus in LAA. Doing so may result in serious patient injury.

Do not use the Clip or applicator to dissect tissue. Doing so may result in tissue damage or tearing. If there are concerns about the presence of adhesions on target anatomy in redo surgery, do not use device. Doing so may result in tissue damage.

⚠️ WARNING ⚠️

This device contains small amounts of Cobalt (CAS# 7440-48-4). Do not use this device if the patient has sensitivity to Cobalt as this may result in an adverse patient reaction.

⚠️ WARNING ⚠️

This device contains nitinol, an alloy of nickel and titanium. Persons with allergic reactions to nickel may suffer an allergic reaction to this implant. Prior to implantation, patients should be counselled on the materials contained in the device, as well as potential for allergy/hypersensitivity.

COMPLICATIONS

Potential complications associated with the use of the AtriClip PRO-Mini LAA Exclusion System and procedure include, but are not limited to, those listed below:

- Air embolism
- Allergic reaction to anesthesia, anticoagulant, implant material
- Anaphylactic shock
- Anesthesia risks
- Aneurysm
- Angina
- Arrhythmia needing medical treatment (new onset)
- Arterial or venous dissection and/or perforation
- Arterial rupture
- Arterial spasm
- Arteriovenous fistula
- Atelectasis (major lung collapse with significant symptoms such as cyanosis, extreme shortness of breath, dyspnea, and/or stabbing pain on the affected side)
- Atrial rupture
- Atrio-esophageal fistula
- AV block requiring permanent pacemaker (new onset)
- Bleeding requiring intervention
- Blood vessel damage
- Cardiac perforation
- Cardiac tamponade
- Cardiac valve injury
- Cerebrovascular accident (CVA)/Transient Ischemic Attack (TIA)/stroke (ischemic or hemorrhagic)
- Chest pain/discomfort
- Compression of coronary artery
- Conduction disturbances
- Congestive heart failure (new onset or exacerbation)
- Coronary artery injury
- Death
- Device breakage/inability to remove
- Device-related death
- Diaphragmatic paralysis (unilateral or bilateral)
- Drug reaction (significant reaction to any procedure related medications requiring treatment, including allergic reaction and anaphylactic shock)
- Emergency during procedure requiring a change in planned access
- Empyema
- Endocarditis (bacterial)
- Esophageal injury
- Esophageal rupture
- Extension of cardiopulmonary/ extracorporeal bypass
- Fever
- Gastric motility disorders
- Gastro-intestinal bleed
- Hematoma
- Hematuria
- Hemothorax
- Hypertension
- Hypotension
- Iatrogenic atrial flutter
- Iatrogenic lung injury (e.g., chest tube placement)
- Ischemia
- Kinking of coronary artery
- LAA dehiscence
- LAA tears
- Left atrial embolism
- Myocardial infarction (MI)
- Nerve injury (phrenic, laryngeal, thoracic, etc.)
- Pain/discomfort
- Pericardial effusion
- Pericarditis
- Permanent pacemaker
- Persistent chest pain (post discharge surgical incision pain, not angina)
- Phrenic nerve paralysis
- Pleural effusion
- Pneumonia
- Pneumothorax
- Postoperative embolic complications
- Pseudoaneurysm
- Pulmonary edema
- Pulmonary embolism
- Renal insufficiency or failure
- Respiratory distress or failure (breathing problems)
- Sepsis
- Stenosis of left circumflex artery
- Sterility-related infection
- Superficial wound infection
- Surgical site infection
- Systemic adverse reaction due to device corrosion
- Thrombus and/or thromboembolism (including deep vein thrombosis)
- Tissue injury
- Tissue perforation
- Tracheal esophageal trauma
- Vascular access complications

INSTRUCTIONS FOR USE

Surgeon judgment, with the assistance of the Guide, should determine what size Clip to apply. This IFU is designed to assist in using this product. It is not a reference to surgical techniques.

CLIP SELECTION

⚠️ WARNING ⚠️

Carefully consider any presurgical treatment the patient may have undergone when selecting Clip size. Preoperative radiotherapy may result in changes to tissue. These changes may, for example, cause the tissue thickness to exceed the indicated range for the selected Clip size. Failure to correctly size the Clip may result in: tissue trauma, dehiscence, tissue tearing, displacement, lack of desired hemostasis, and/or incomplete exclusion of the structure.

- Using the Guide, determine correct selection of the AtriClip Mini (See Figure 2). Clip sizes are located on the device package.

Labeled Clip Size	LAA Size Range
35 mm	29 – 35 mm (1.14 – 1.38 in)
40 mm	34 – 40 mm (1.34 – 1.57 in)
45 mm	39 – 45 mm (1.54 – 1.77 in)
50 mm	44 – 50 mm (1.73 – 1.97 in)

⚠️ WARNINGS ⚠️

Do not use on a compressed LAA less than 29 mm (1.14 in) in width and 1 mm (0.04 in) wall thickness. Doing so may result in: tissue trauma, dehiscence, tissue tearing, displacement, and/or lack of desired hemostasis. Do not use on a compressed LAA greater than 50 mm (1.97 in). Doing so may result in incomplete exclusion of the structure.

- Using sterile technique, remove the AtriClip PRO-Mini LAA Exclusion System from its packaging.

⚠️ WARNING ⚠️

If the sterile package is damaged and/or the sterile barrier is breached, discard device and DO NOT USE to avoid the risk of patient infection.

⚠️ **CAUTION:** Do not drop the device as this may induce damage to the device. If the device is dropped, do not use. Replace with a new device.

- Using the Activation Lever on the handle, gently open and close the Clip to assure proper function.

ARTICULATION OF END EFFECTOR

4. By pushing down and pulling the Articulation Lock backwards (proximal), the End Effector of the AtriClip PRO-Mini LAA Exclusion System may be manually articulated up and down and side-to-side by either the articulation levers or pressing on the End Effector. The Clip and End Effector can articulate 30° left or right and 30° up or down to take into account anatomical variations in the patient's anatomy (See Figure 3).

5. To lock the End-Effector in position, press the Articulation Lock forward.

⚠️ **CAUTION:** Do not attempt to articulate the End Effector while in the locked position. Force applied while in the locked position may cause damage to the device.

CLIP POSITIONING

⚠️ WARNINGS ⚠️

Take care not to damage surrounding structures when maneuvering and positioning the Clip. Position and deploy Clip in a manner that provides direct visualization of all tissues being accessed. Direct visualization, in this context, requires that the surgeon is able to see the heart directly, with or without assistance from a camera, endoscope, etc., or other appropriate viewing technologies. Poor visualization may result in suboptimal placement and damage or obstruction of surrounding structures.

⚠️ **CAUTION:** Do not kink or bend the Shaft as this may affect device performance.

- Maneuver the AtriClip PRO-Mini LAA Exclusion System into the targeted dissection plane.
- Gently open the Clip by squeezing the Activation Lever.

NOTE: The Clip automatically locks in the fully open position by means of a Locking Trigger on the handle of the device. The lock can be disengaged by pressing the Lever Release Trigger.

- Orient the Clip applicator with preloaded Clip at the tip of the LAA (See Figure 4).
- Gently position the Clip at the base of the LAA (See Figure 5).
- Position the Clip in a manner that provides clear visualization of all tissues being accessed.
- While the Clip is still affixed to the Deployment Device, ensure that no surrounding structures interfere with or are damaged by the Clip, and that the Clip is placed correctly.
- If the Clip is not placed correctly, gently open the Clip and reposition as needed.
- After the Clip is positioned correctly, grasp the Activation Lever, depress the Lever Release Trigger, and slowly release the Activation Lever to allow the Clip to close.

DEPLOYMENT

⚠️ WARNING ⚠️

Carefully evaluate Clip position, tissue thickness, and tissue width prior to Clip deployment. To determine appropriate Clip size, refer to the Guide Instructions for Use. Failure to correctly size or deploy the Clip may result in: tissue trauma, dehiscence, tissue tearing, displacement, and/or lack of desired hemostasis. Unless medically necessary, do not attempt to reposition or remove the Clip after deployment. This may result in tissue damage or tearing.

- Deploy the Clip by slowly pulling the Deployment Tab at the proximal end of the Handle.

NOTE: The Deployment Tab with steel cables may be completely removed from the end of the Handle (See Figure 6).

⚠️ **CAUTION:** Take care to minimize manipulation of the LAA and Clip after Clip deployment.

- Unlock the End Effector articulation and carefully remove it from the LAA, leaving the Clip and attachment suture behind (See Figure 7).

NOTE: After pulling the Deployment Tab, the AtriClip PRO-Mini LAA Exclusion System cannot be used to reposition or remove the Clip.

DISPOSAL INFORMATION

After use, this device should be treated as medical waste and disposed of following hospital protocol.

SSCP

A summary of the safety and clinical performance of the device can be found in the European database on medical devices (Eudamed) at <https://ec.europa.eu/tools/eudamed> by entering the Basic UDI-DI associated with the device.

Product Code(s)	Basic UDI-DI
PROM35, PROM40, PROM45, PROM50	084014390000000000016ZQ

SERIOUS INCIDENT

Any serious incident that has occurred in relation to this device should be reported to AtriCure and the competent authority of the Member State in which the user and/or patient is located.

IFU-0551.A
2026-06
Page 2 of 8

RETURN OF USED PRODUCT

If, for any reason, this product must be returned to AtriCure, Inc., a return goods authorization (RGA) number is required from AtriCure, Inc., prior to shipping.

If the product has been in contact with blood or body fluids, it must be thoroughly cleaned and disinfected before packing. It should be shipped in either the original carton or an equivalent carton, to prevent damage during shipment; and it should be properly labeled with an RGA number and an indication of the biohazardous nature of the contents of shipment.

Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number may be obtained from AtriCure, Inc.

DISCLAIMER STATEMENTS

Users assume responsibility for approving the acceptable condition of this product before it is used, and for ensuring that the product is only used in the manner described in these instructions for use, including, but not limited to, ensuring that the product is not re-used.

Under no circumstances will AtriCure, Inc. be responsible for any incidental, special or consequential loss, damage, or expense, which is the result of the deliberate misuse or re-use of this product, including any loss, damage, or expense which is related to personal injury or damage to property.

Not all products/models may be available in your specific country. Contact your local AtriCure representative to check availability.

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MRI SAFETY INFORMATION

⚠ MR CONDITIONAL

The AtriClip Mini is MR Conditional. A patient with the AtriClip Mini may be safely scanned under the following conditions. Failure to follow these conditions may result in injury to the patient.

- Static Magnetic Field Strength (T):** 1.5-Tesla or 3.0-Tesla
- Static Magnetic Field Orientation:** Horizontal
- Maximum Spatial Field Gradient (T/m and gauss/cm):** 40-T/m (4,000-gauss/cm)
- RF Excitation Polarization:** Circularly Polarized (CP) (i.e., Quadrature-Transmission)
- Transmit RF Coil:** Any transmit RF coil may be used
- Receive RF Coil:** Any receive-only RF coil may be used
- MR System Operating Mode:** Normal Operating Mode
- Maximum Whole Body Averaged SAR (W/kg):** 2-W/kg (Normal Operating Mode)
- Scan Duration and Wait Time:** Whole body averaged SAR of 2-W/kg for 60 minutes of continuous RF exposure (i.e., per pulse sequence or back-to-back sequences/series without breaks)
- MR Image Artifact:** The presence of the AtriClip Mini may produce an MR image artifact. Imaging protocol modifications may be necessary to compensate for the MR image artifact.

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AtriCure Inc.
AtriCure Europe B.V.