

**EU Indications:** AtriCure's cryoICE cryoSPHERE cryoablation probes are sterile, single use devices intended for use in blocking pain by temporarily ablating intercostal peripheral nerves.

Please review the Instructions for Use for a complete listing of contraindications, warnings, precautions and potential adverse events prior to using these devices.

*cryo***SPHERE**<sup>®</sup>

The Coolest Innovation  
in Pain Management.

## *cryo***SPHERE**<sup>®</sup>

- Enhanced Flexibility
- Uniform Application
- Engineered for Intercostal Nerve Cryoablation

### **Indications**

AtriCure's cryoICE<sup>®</sup> cryoSPHERE<sup>®</sup> cryoablation probes are sterile, single use devices intended for use in blocking pain by temporarily ablating intercostal peripheral nerves.

### **ATRICURE, INC.**

7555 Innovation Way  
Mason, Ohio 45040 USA  
Customer Service: +1 (866) 349-2342  
CustomerService@AtriCure.com  
[www.AtriCure.com](http://www.AtriCure.com)

### **ATRICURE EUROPE B.V.**

De Entree 260  
1101 EE Amsterdam  
Netherlands  
Tel: +31 (0) 20-7005560  
Fax: +31 (0) 20-7005561  
SalesSupportEU@AtriCure.com  
[www.AtriCure.com/International](http://www.AtriCure.com/International)


PM-EU-1697B-0524-G

**AtriCure**

**AtriCure**



# The Coolest Innovation in Pain Management.



## Features

**8mm Ball Tip** Designed for minimally invasive applications and placement on intercostal nerve

**Bendable Distal Shaft** Bendable distal shaft up to 180 degrees


**Braided Insulation** Technology that is designed to maintain a temperature warmer than -10°C in tissue

**Reinforced Proximal Shaft** Provides rigidity at the level of the incision points

**Ergonomic Handle Design** Provides ease of use and comfort during cryoablation

**Active Defrost** Allows probe to be easily removed from cryoablation site

## Malleable Cryoablation Probe



## cryoICE BOX® - V6

ACM1 (US)



## cryoICE Cart & Upgrade Kits

- ASCS (cart with double tank holder, US only)
- ASCC (single tank holder upgrade)
- ASCC2 (double tank holder upgrade)

20 lbs N2O tanks not provided by AtriCure



## Foot Switch (optional)

For use with ACM1/ACM2  
FTSWITCH

