

# AtriClip PRO2 Device

AtriClip PRO Devices can be used in totally thoracoscopic procedures in conjunction with other heart heart surgery.



## Features

- Hoopless End Effector
- Small Diameter End Effector (12 mm)
- Active Articulation Levers

LAA Exclusion System		
Size	AtriClip PRO Device	AtriClip PRO2 Device
35 mm	PRO135	PRO235
40 mm	PRO140	PRO240
45 mm	PRO145	PRO245
50 mm	PRO150	PRO250

The AtriClip LAA Exclusion System is indicated for the occlusion of the heart's left atrial appendage.

Please review the Instructions for Use for a complete listing of contraindications, warnings, precautions and potential adverse events prior to using these devices.

# AtriClip PRO2<sup>®</sup> Device

## A Simplified Approach to LAA Exclusion



**AIRICURE EUROPE B.V.**  
De Entree 260  
1101 EE Amsterdam  
Pays-Bas  
Tél. : +31 (0) 20-7005560  
Fax : +31 (0) 20-7005561  
SalesSupportEU@AtriCure.com  
europe.atricure.com  
PM-EU-2371A-1224-G

**AtriCure**

**AtriCure**


# AtriClip PRO2 Device

**AtriClip PRO<sup>®</sup>**  
Devices can now  
be used in totally  
thoracoscopic  
procedures.



Greater than 30°  
Pitch and Yaw

Rigid Shaft



12 mm  
Diameter  
End  
Effector

