

# AtriClip FLEX•V® Device

# AtriCure


180° of Shaft  
Rotation



New Clip Opening  
Handle with Reduced  
Fatigue Design



Malleable Shaft



Enhanced LAA Visualization

One-Handed Application • Suture-Less Clip Deployment

Small Footprint to Minimize the Interference with Adjacent Structures

The AtriClip FLEX•V Device Applies Continual Force  
to the Base of the LAA as it Atrophies


Open-Ended  
AtriClip

Tip-First  
Closure

Continual Force  
Applied to the  
Base of the LAA



AtriClip® Device Performance  
with a Clip Deployment Trigger



The AtriClip™ LAA Exclusion System is indicated for the occlusion of the heart's left atrial appendage.

Please review the Instructions for Use for a complete listing of contraindications, warnings, precautions and potential adverse events prior to using these devices.

**ATRICURE EUROPE B.V.**

De Entree 260 1101 EE Amsterdam Pays-Bas

Tél. : +31 (0) 20-7005560

Fax : +31 (0) 20-7005561

SalesSupportEU@AtriCure.com / europe.atricure.com

PM-EU-2374A-1224-G

**AtriCure**