

SEIZE THE CHANCE, EXCLUDE THE LAA

AtriCure's AtriClip products are the most widely sold Left Atrial Appendage (LAA) management devices worldwide.

>97%
successful LAA
exclusion^{1*}

Epicardial Exclusion

- Implant is not in the blood stream
- Ischemic injury electrically isolates the LAA
- Atrophy and resorption of the LAA

Dynamic Closing Force

Continuous closing force maintains LAA exclusion throughout changes to the tissue caused by ischemia

Parallel / Linear Closing

Minimized occurrence of tissue folds with optimal apposition of tissue along long axis of LAA ostia

Tissue Compression / Atraumatic

No cutting, non-piercing, and atraumatic compression reduces risk of tissue tearing and bleeding



AtriCure


The closed loop clips

Seize the chance, exclude the LAA

AtriClip Device Systems | Open accesses

ATRICLIP FLEX

Features




The open clips

Seize the chance, exclude the LAA

AtriClip Device Systems | Open accesses

ATRICLIP FLEX•V


Features



AtriClip Device Systems | Minimal accesses

ATRICLIP PRO

Features



AtriClip Device Systems | Minimal accesses

ATRICLIP PRO•V

Features



ATRICLIP PRO 2

Features



AtriClip FLEX

Device	Size
ACH235	35 mm
ACH240	40 mm
ACH245	45 mm
ACH250	50 mm

AtriClip FLEX•V

Device	Size
ACHV35	35 mm
ACHV40	40 mm
ACHV45	45 mm
ACHV50	50 mm

AtriClip PRO

Device	Size
PRO135	35 mm
PRO140	40 mm
PRO145	45 mm
PRO150	50 mm

AtriClip PRO•V

Device	Size
PROV35	35 mm
PROV40	40 mm
PROV45	45 mm
PROV50	50 mm

Selection device
CGC

AtriClip PRO 2

Device	Size
PRO235	35 mm
PRO240	40 mm
PRO245	45 mm
PRO250	50 mm

The AtriClip LAA Exclusion System is indicated for the occlusion of the heart's left atrial appendage.

1. Benussi, S. et al, Thoracoscopic Appendage Exclusion With an AtriClip Device As a Solo Treatment for Focal Atrial Tachycardia. Circulation, 2011, 123:1575-8.

ATRICURE EUROPE B.V.

De Entree 260
1101 EE Amsterdam
Netherlands
Tel: +31 (0) 20-7005560
Fax: +31 (0) 20-7005561
SalesSupportEU@AtriCure.com
europe.atricure.com

© AtriCure, Inc.
PM-EU-3094A-1124-G

AtriCure