



A comprehensive
surgical AF strategy
combining a box lesion
and LAA exclusion.

Product	Code	Information
Isolator Synergy EnCompass Clamp	OSH	Includes GPM100
AtriClip FLEX-Mini / FLEX•V	ACHM35 / ACHV35	35 mm
AtriClip FLEX-Mini / FLEX•V	ACHM40 / ACHV40	40 mm
AtriClip FLEX-Mini / FLEX•V	ACHM45 / ACHV45	45 mm
AtriClip FLEX-Mini / FLEX•V	ACHM50 / ACHV50	50 mm
AtriClip Selection Guide (Sizer)	CGG	
Glidepath Magnetic Guide (Clamp)	GPM100	Individually package

The AtriClip LAA Exclusion System is indicated for use in patients at high risk of thromboembolism for whom left atrial appendage exclusion is warranted.

The Isolator Synergy EnCompass Clamp is indicated for the ablation of cardiac tissue for treatment of cardiac arrhythmias, including atrial fibrillation.

Please review the Instructions for Use for a complete listing of contraindications, warnings, precautions and potential adverse events prior to using these devices.

*Data on file. Ablations performed on ex-vivo bovine myocardium and epicardium, where transmuralit is defined as lesion depth equals tissue thickness. The transmuralit tests were performed at 2, 5, 10 and 15mm of tissue thickness.

¹Mean closure rate among the 3 studies is 98.1% per data on file at AtriCure Inc. Ailawadi G. et al. J Thorac Cardiovasc Surg 2011;142:1002–9, 1009.e1; Gerdisch MW. et al. Innovations (Phila). 2022;17(6):463-470; Cartledge R et al. Interactive CardioVascular and Thoracic Surgery 34 (2022) 548–555.

²Starck, C.T. et al. (2012). Epicardial left atrial appendage clip occlusion also provides the electrical isolation of the left atrial appendage. Interactive cardiovascular and thoracic surgery, 15(3):416-18. This human clinical study demonstrated acute electrical isolation. Chronic electrical isolation has not been evaluated in human clinical studies.

ATRICURE EUROPE B.V.
De Entree 260 1101 EE Amsterdam Pays-Bas
Tél. : +31 (0) 20-7005560
Fax : +31 (0) 20-7005561
SalesSupportEU@AtriCure.com
europe.atricure.com

PM-EU-3687B-0928-G

AtriCure



BoxX Procedure

Ablations Simplified, Results Amplified

Transmural Ablation

Designed to create contiguous, transmural ablations across varying tissue thicknesses.

Intelligent Energy Delivery


The Synergy algorithm measures tissue response 50 times per second adjusting energy delivery throughout each ablation.

Parallel Closure & Uniform Pressure

Designed to provide consistent tissue contact along the ablation line.

Guided Placement

The Glidepath™ Magnetic Guide supports controlled positioning of the Isolator Synergy EnCompass® clamp around the target anatomy.




Isolator Synergy EnCompass Clamp

The EnCompass clamp is designed to make surgical ablation transmural, contiguous and easy to perform for non-atriotomy patients.



EnCompass Clamp gives you:

- Posterior wall access
- Procedural flexibility
- Minimal tissue dissection




The Glidepath Magnetic Guide, designed to safely place the clamp



AtriClip® Devices


Proven Performance. Expanded Possibilities:
Two distinct designs to support different patients, procedures and surgical preferences.

AtriClip Flex-Mini® Device



Small Footprint to Minimize the Interference with Adjacent Structures and Enhance LAA Visualization

AtriClip FLEX-V® Device



Seize the chance, exclude the LAA



Stroke Risk Reduction*

LAAOS III demonstrated a 33% reduction in stroke or systemic embolism with surgical LAA occlusion, with AtriClip device among the techniques used. Long-term AtriClip device data further support durable LAA exclusion.

Epicardial Exclusion

Excludes the LAA without leaving an implant in the bloodstream, while ischemic injury provides electrical isolation of the LAA.²

Dynamic Closing Force

Continuous closing force maintains LAA exclusion as tissue changes.

Atraumatic Tissue Compression

No cutting or piercing of the LAA, reducing the risk of tissue tearing and bleeding.

*Originated from the LAA