

AtriClip PRO® Device

Complete LAA Occlusion

>98%
successful LAA
exclusion¹

AtriCure's AtriClip products are the most widely sold Left Atrial Appendage (LAA) management devices worldwide.

Epicardial Exclusion

- Implant is not in the blood stream
- Ischemic injury electrically isolates the LAA
- Atrophy and resorption of the LAA

Dynamic Closing Force

Continuous closing force maintains LAA exclusion throughout changes to the tissue caused by ischemia

Parallel / Linear Closing

Minimized occurrence of tissue folds with optimal apposition of tissue along long axis of LAA ostia

Tissue Compression / Atraumatic

No cutting, non-piercing, and atraumatic compression reduces risk of tissue tearing and bleeding

Stroke* risk reduction

Indicated for use in patients at high risk of thromboembolism for whom left atrial appendage exclusion is warranted

*Originated from the LAA



Now recommended **Class I**
by ESC/EACTS 2024
Recommendations

AtriClip PRO Devices can be used in totally thoracoscopic procedures in conjunction with other heart surgery.

Minimal accesses

AtriClip PRO

Features

- $\pm 30^\circ$ omni directional head articulation with lock
- Quick deploy feature
- 25 cm Rigid Shaft
- Articulation Lock
- Quick Release
- Clip Opening Lever



Device	Size
PRO135	35 mm
PRO140	40 mm
PRO145	45 mm
PRO150	50 mm

Selection device
CGG

The AtriClip LAA Exclusion System is indicated for use in patients at high risk of thromboembolism for whom left atrial appendage exclusion is warranted.

Please review the Instructions for Use for a complete listing of contraindications, warnings, precautions and potential adverse events prior to using these devices.

1. Mean closure rate among the 3 studies is 98.1% per data on file at AtriCure, Inc. Ailawadi G. et al. J Thorac Cardiovasc Surg 2011;142:1002-9, 1009.e1; Gerdisch MW. et al. Innovations (Phila). 2022;17(6):463-470; Cartledge R et al., Interactive CardioVascular and Thoracic Surgery 34 (2022) 548-555

ATRICURE EUROPE B.V.

De Entree 260
1101 EE Amsterdam
Netherlands
Tel: +31 (0) 20-7005560
Fax: +31 (0) 20-7005561
SalesSupportEU@AtriCure.com
Europe.AtriCure.com

© AtriCure, Inc.
PM-EU-3955A-0127-G

AtriCure