

AtriClip PRO-Mini® Device

Optimized visualization and placement
for minimally invasive procedures

>98%
successful LAA
exclusion¹

AtriCure's AtriClip products are the most widely sold Left Atrial Appendage (LAA) management devices worldwide.

Epicardial Exclusion

- Implant is not in the blood stream
- Ischemic injury electrically isolates the LAA
- Atrophy and resorption of the LAA

Dynamic Closing Force

Continuous closing force maintains LAA exclusion throughout changes to the tissue caused by ischemia

Parallel / Linear Closing

Minimized occurrence of tissue folds with optimal apposition of tissue along long axis of LAA ostia

Tissue Compression / Atraumatic

No cutting, non-piercing, and atraumatic compression reduces risk of tissue tearing and bleeding

Stroke* risk reduction

Indicated for use in patients at high risk of thromboembolism for whom left atrial appendage exclusion is warranted



Class I Recommendation
ESC/EACTS 2024 Guidelines

*Originated from the LAA

Minimal accesses

AtriClip PRO-Mini

AtriClip Mini for Minimally Invasive
Surgical Procedures:

- Clinically proven to provide both mechanical and electrical isolation of the LAA^{2,3}
- Durable and atraumatic LAA exclusion^{2,4}
- Same fabric as existing AtriClips with a tighter weave that is form fitting in nature for optimized visualization and a smaller footprint
- Custom fit to specific anatomy with multiple sizes available
- 12 mm port compatibility


Combined with Proven PRO Features:

- Hoopless End Effector
- Active Articulation Levers and Locking Mechanism
- Quick Release
- Ambidextrous Lever Release Trigger
- Clip Opening Lever
- Rigid Shaft Greater than 30°Pitch and Yaw
 - Right - Up and Down (Pitch)
 - Left - Side to Side (Yaw)

Constant innovation resulting
in three generations of AtriClip



AtriClip PRO-Mini is preloaded
with the AtriClip Mini that is
60% smaller profile*.



*than the next lowest profile on the market.

AtriClip PRO-Mini

Device	Size
PROM35	35 mm
PROM40	40 mm
PROM45	45 mm
PROM50	50 mm

Selection device
CGG

The AtriClip LAA Exclusion System is indicated for use in patients at high risk of thromboembolism for whom left atrial appendage exclusion is warranted.

Please review the Instructions for Use for a complete listing of contraindications, warnings, precautions and potential adverse events prior to using these devices.

1. Mean closure rate among the 3 studies is 98.1% per data on file at AtriCure, Inc. Ailawadi G. et al. J Thorac Cardiovasc Surg 2011;142:1002–9, 1009.e1; Gerdisch MW. et al. Innovations (Phila). 2022;17(6):463-470; Cartledge R et al., Interactive CardioVascular and Thoracic Surgery 34 (2022) 548–555
2. Data on file.
3. Starck, C.T., Steffel, J., Emmert, M.Y. et al. (2012). Epicardial left atrial appendage clip occlusion also provides the electrical isolation of the left atrial appendage. Interactive cardiovascular and thoracic surgery, 15(3):416-418. This human clinical study demonstrated acute electrical isolation. Chronic electrical isolation has not been evaluated in human clinical studies.
4. Emmert, M.Y., Puippe, G., Baumüller, S. et al. (2014). Safe, effective and durable epicardial left atrial appendage clip occlusion in patients with atrial fibrillation undergoing cardiac surgery: first long-term results from a prospective device trial. Eur J Cardiothorac Surg, 45(1):126-131.

ATRICURE EUROPE B.V.

De Entree 260
1101 EE Amsterdam
Netherlands
Tel: +31 (0) 20-7005560
Fax: +31 (0) 20-7005561
SalesSupportEU@AtriCure.com
Europe.AtriCure.com

© AtriCure, Inc.
PM-EU-4502B-0728-G

AtriCure